Histological Clearing Agent

Xylene Substitute

(Issue Date: 06/24/2022)

Clearify™, The product and its purpose:

Clearify™ is intended to be used as a xylene substitute in tissue processing and staining of histological and cytological slides preparations. **Clearify™** can be used both in the clearing and deparaffinization steps for histological and clearing steps for cytological slide preparations allowing for further qualitative in-vitro diagnostic examinations in the field of cellular pathology. **Clearify™** may also be used as a clearing solvent in both the manual and automated tissue processing procedures. This device is intended to be used by trained laboratory professionals.

Clearify™ Ready to Use

Principles and Procedures:

Before the staining process, tissue samples must be prepared through the following stages: Fixation, dehydration, embedding, and sectioning. Fixation preserves the structure of the tissue in its natural form. During fixation with the use of Neutral Buffered Formalin, **Clearify™** is used after the last alcohol stage and before immersing the tissue specimen in liquid paraffin. Dehydration is the process of removing water from the sample by adding Ethanol. After completion of dehydration, Xylene is used to remove the Ethanol. Embedding is the process of putting the sample into a paraffin wax to enhance the process of extracting cellular structures. **Clearify™** is used to remove the paraffin from thin paraffin sections before staining and to prepare the sections for the descending alcohol series. Sectioning involves mounting the specimen on a microtome and cutting it into sections in order to stain the specimen for examination.

Histoprocessing:

Neutral Buffered Formalin
Ethanol 100% Denatured
Clearify™
Paraffin at 60 °C

Clearify™, Instructions for use:

Clearify™ requires small adjustments to routine procedures. Rotate and change the final tissue processor alcohol frequently to ensure no water is carried over. Two changes of **Clearify™** at one hour each are needed to insure complete clearing. For best results use **Clearify™** with Alcohols denatured with Isopropyl Alcohol, but not Methanol. **Clearify™** is not recommended for paraffin cleaning cycles on tissue processors.

Deparaffinizing:

1. Remove slides from oven and place into **Clearify™** while paraffin is molten.
2. Deparaffinize through two changes of **Clearify™** for 10 minutes each.

Dehydrating, Clearing, and Coverslipping

1. Dehydrate slide through 4 changes of Absolute Alcohol.
2. Clear slides through 3 changes of **Clearify™**. 2 minutes each in the first two changes and 5 minutes in the final change.

Clearify™, Storage and stability:

Clearify™ should be stored at 15-30 °C and is stable under normal laboratory conditions. See product label for shelf life. Discard product if leaking or if solution appears cloudy discontinue use.

Clearify™ Shelf Life:

After first opening, **Clearify™** can be used until the stated expiry date when stored. Keep container tightly closed in 15-30 °C temperature conditions.

Clearify™, Precautions:

See Safety Data Sheet (SDS) for warnings and precautions.

Clearify™ Main Product Components:

Naphtha Isoparaffinic Hydrocarbon
CAS No.: 64742-48-9
EC No.: 265-150-3

CLEARIFY™

Histological Clearing Agent Xylene Substitute

Clearify™, Disposal:

Dispose **Clearify™** in accordance with all local, regional, national, territorial, provincial, international regulations. Within the EU the currently applicable REGULATION (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006 applies.

Hazard Classification

Cat. No. CACLELT, CACLEGAL

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and / or its authorized representative and to your national authority.

Refer to product label and information stated in the Safety Data Sheet.

The safety data sheet is available upon request and on www.statlab.com

Literature:

Histotechnology 3rd Edition, 2009

Freida L Carson, Christa Haldik

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	SLMP LLC d.b.a. StatLab Medical Products 2090 Commerce Drive McKinney, Texas 75069 USA Tel: 1-800-442-3573
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	 GHS08	Danger - Health Hazard	

**Copyright © 2019-2021 SLMP, LLC, DBA
StatLab. All rights reserved.**

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.

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Copyright © 2019-2021 StatLab Medical Products.
Copyright © 2019-2021 StatLab Medical Products. Tous droits réservés.
Copyright © 2019-2021 StatLab Medical Products. Med ensamrätt.
Copyright © 2019-2021 SatLab Medical Products. Всички права запазени.
Copyright © 2019-2021 StatLab Medical-produkter. Alle rettigheder forbeholdes.
Copyright © 2019-2021 StatLab Medical Products. Všechna práva vyhrazena.