American MasterTech
scientific laboratory supplies

ExCell Plus

Ready to Use Histological Tissue Fixative

(Rev. 3/8/18)

ExCell Plus™, The product and its purpose:

ExCell Plus™ is a Histological tissue fixative formulated as a replacement for 10% Formalin and intended for In Vitro Diagnostic use. Tissue fixed with **ExCell Plus™** can be used for routine Special Stains and Immunohistochemical stains (IHC). Antigen retrieval will likely still be needed for IHC staining. **ExCell Plus™** is compatible with other routine fixatives and alcohols. Active ingredient include; Glyoxal (10-20%) and Ethanol (1-5%).

ExCell Plus™, Instructions for use:

Use **ExCell Plus™** as you would 10% Formalin. Most surgical tissue specimens require 6 to 24 hours for complete fixation. Large tissues, whole organs, mastectomies, and placentas require at least 72 hours (bisecting will help ensure complete penetration and fixation).

ExCell Plus™, Tissue storage:

Tissue specimens may be stored in **ExCell Plus™** for several months without tissue degradation or adverse staining. **ExCell Plus™** will not over-harden tissue specimens, even when stored for long periods of time.

ExCell Plus™, Storage and Stability:

ExCell Plus™ should be stored at 15-30 °C and is stable under normal laboratory conditions. See product label for shelf life. Color of **ExCell Plus™** may vary slightly but this does not affect product performance or stability. If **ExCell Plus™** freezes during transit, let product thaw at room temperature until completely liquefied; performance and stability of **ExCell Plus™** is unaffected.

ExCell Plus™, Precautions:

See Safety Data Sheet (SDS) for warnings and precautions.

ExCell Plus™, Disposal:

Dispose **ExCell Plus™** in accordance with all local, state, federal, environmental and/or EU regulations.

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CE MARKINGS AND DESIGNATIONS:

	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
	Batch Code		Use By		Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use			
	Irritant		Health Hazard			

Copyright © 2015-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.