Human Colon Tumor

## Alcian Blue pH 1.0 Stain Kit Procedure

**100ml Kit Item #:** KTABP1

**Liter Kit Item#:** N/A

**Pint Kit Item #:** KTABP1PT

**Gallon Kit Item#:** N/A

**Control Slide(s)**
**Item#**
**Included Components**

Colon

CSC1025P

3% Acetic Acid

Small Intestine

CSS0325P

Alcian Blue Stain pH 1.0

Nuclear Fast Red Stain

**PRINCIPLE AND RESULTS:** This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify strongly sulfated mucosubstances. Strongly sulfated mucosubstances are stained blue, cytoplasm and background pale pink, and nuclei pink to red.

**SPECIMEN CRITERIA:** Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

**STORAGE AND USAGE NOTES:** Store/Use each component according to the temperature and expiration on the label.

**PRECAUTIONS:** For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

### STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

| # | Action        | With                            | Heat °C | Time |      | Details                                                                   |
|---|---------------|---------------------------------|---------|------|------|---------------------------------------------------------------------------|
|   |               |                                 |         | Mins | Secs |                                                                           |
| 1 | Deparaffinize | Xylene or Substitute, 2 changes | --      | 5    | --   | 5 minutes each change or as required if using a xylene substitute.        |
| 2 | Rinse         | Absolute Alcohol, 3 changes     | --      | 1    | --   | 1 minute each change or as required if using graded alcohols.             |
| 3 | Rinse         | Running Tap Water               | --      | 1    | --   |                                                                           |
| 4 | Immerse       | 3% Acetic Acid                  | --      | 3    | --   | Without rinsing, continue to next step.                                   |
| 5 | Immerse       | Alcian Blue Stain pH 1.0        | --      | 30   | --   | Once complete, <b>rinse in running tap water (1 minute)</b> and continue. |
| 6 | Immerse       | Nuclear Fast Red Stain          | --      | 5    | --   | Once complete, <b>rinse in running tap water (1 minute)</b> and continue. |
| 7 | Dehydrate     | Absolute Alcohol, 3 changes     | --      | 1    | --   | 1 minute each change.                                                     |
| 8 | Clear         | Xylene or Substitute, 3 changes | --      | 1    | --   | 1 minute each change or as required if using a xylene substitute.         |
| 9 | Coverslip     | Permanent Mounting Media        | --      | --   | --   |                                                                           |

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology, 1980, 172 – 173
2. With modifications by AMTS R&D Department, 1979-2019.

### AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. A minimum of 4 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 15 baths to run the complete procedure.

**TEST YIELD:** \*Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

| Bath Type              | Uses | Slides | S.C. | Bath Type              | Uses | Slides | S.C. |
|------------------------|------|--------|------|------------------------|------|--------|------|
| 20ml Plastic Slide Jar | 25   | 500    | 20   | 250ml Glass Stain Dish | 2    | 375    | 188  |
| 30ml Glass Coplin Jar  | 16   | 480    | 30   | 200ml Bath Autostainer | 2    | 320    | 160  |
| 40ml Hellendahl Jar    | 12   | 480    | 40   | 400ml Bath Autostainer | 1    | 340    | 340  |

### CE MARKINGS AND DESIGNATIONS:

|            |                                    |                                                                                   |                                   |                                                                                            |                |                                                                                                                  |
|------------|------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------|--------------------------------------------------------------------------------------------|----------------|------------------------------------------------------------------------------------------------------------------|
| <b>REF</b> | Catalogue Number                   |  | Temperature Limitation            |           | Manufacturer   | American MasterTech Scientific<br>1330 Thurman St.<br>Lodi, CA 95240 USA<br>Tel 800 860 4073<br>Fax 209 368 4136 |
| <b>LOT</b> | Batch Code                         |  | Use By                            | <b>EC REP</b>                                                                              | Representative | Emergo Europe<br>Prinsessegracht 20<br>2514 AP The Hague<br>The Netherlands                                      |
| <b>IVD</b> | In Vitro Diagnostic Medical Device |  | Consult Instructions Prior to Use | <br>GHS08 | Health Hazard  |                               |

### CONTACT INFORMATION:

**American MasterTech Scientific**  
 1330 Thurman St.  
 Lodi, CA 95240 USA  
 Tel 800 860 4073  
 Fax 209 368 4136  
[www.americanmastertech.com](http://www.americanmastertech.com)

**StatLab**  
 2090 Commerce Drive  
 McKinney, TX 75069 USA  
 Tel 800 442 3573  
 Fax 972 436 1369  
[www.statlab.com](http://www.statlab.com)

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology, 1980, 172 – 173
2. With modifications by AMTS R&D Department, 1979-2019.

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**COMPOSANTS INCLUS:** 3% Acetic Acid, Alcian Blue Stain pH 1.0, Nuclear Fast Red Stain

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