Skeletal Muscle

## Myofibrillar ATPase Stain Kit Procedure

100ml Kit Item #: KTATP

Liter Kit Item#: KTATPLT

Pint Kit Item #: KTATPPT

Gallon Kit Item#: N/A

Control Slide(s) | Item#

N/A

### Included Components

Type 1 Buffer	Cobalt Chloride Solution
Type 1&2BC Buffer	1N Sodium Hydroxide
Alkaline Buffer Stock	1N Hydrochloric Acid
Calcium Chloride Solution	ATP Powder

### Components Not Included

Ammonium Sulfide (Item#: SPA2226)

**PRINCIPLE AND RESULTS:** This kit is intended for use by laboratory professionals to stain routinely prepared frozen sections (in vitro) to identify and differentiate muscle types. Type 1 Buffer stains small type 1 fibers black and type 2C fibers grey. Type 1&2BC Buffer stains small type 1 fibers black, 2B and 2C fibers grey to black. Type 2AB Buffer stains type 2A and 2B fibers black.

**SPECIMEN CRITERIA:** A Minimum of three, Fresh, unfixed, properly prepared and frozen, 12 µm serialized muscle tissue sections. After preparing buffers and immediately before beginning staining procedure, allow frozen sections to warm up to room temperature for up to a maximum of five minutes.

**PRECAUTIONS:** For use by laboratory professionals. Sensitive and properly calibrated pH Meter is required for buffer and working solutions preparations. Dilute Sodium Hydroxide or Hydrochloric Acid pH adjustment solutions with distilled or DI water as necessary to achieve proper pH of working solutions. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

**TYPE 1 BUFFER PREPARATION:** Prepare solution at time of use. Solution expires after 8 hours or out of pH range.

#	Action	Amount	Chemical/Reagent	Details
1	Add	20ml	Type 1 Buffer	Into a chemically clean container or new/unused plasticware.
2	pH	Drops	1N Hydrochloric Acid (to lower pH) or 1N Sodium Hydroxide (to raise pH)	pH must be between 4.25 and 4.35 (4.30 is ideal).

**TYPE 1&2BC BUFFER PREPARATION:** Prepare solution at time of use. Solution expires after 8 hours or out of pH range.

#	Action	Amount	Chemical/Reagent	Details
1	Add	20ml	Type 1&2BC Buffer	Into a chemically clean container or new/unused plasticware.
2	pH	Drops	1N Hydrochloric Acid (to lower pH) or 1N Sodium Hydroxide (to raise pH)	pH must be between 4.55 and 4.65 (4.60 is ideal).

**TYPE 2AB BUFFER PREPARATION:** Prepare solution at time of use. Solution expires after 90 minutes, two uses, or out of pH range.

#	Action	Amount	Chemical/Reagent	Details
1	Add	20ml	Alkaline Buffer Stock	Into a chemically clean container or new/unused plasticware.
2	pH	Drops	1N Hydrochloric Acid (to lower pH) or 1N Sodium Hydroxide (to raise pH)	Starting pH must be 10.45 to allow for natural pH drop during procedure. If two back-to-back uses, re-pH to 10.45 before second batch of slides.

**ATPase INCUBATION SOLUTION PREPARATION:** Prepare solution at time of use. Solution expires after 90 minutes, two uses, or out of pH range.

#	Action	Amount	Chemical/Reagent	Details
1	Add	20ml	Alkaline Buffer Stock	Into a chemically clean container or new/unused plasticware.
2	Add	32mg	ATP Powder	Mix thoroughly.
3	pH	Drops	1N Hydrochloric Acid (to lower pH) or 1N Sodium Hydroxide (to raise pH)	Starting pH must be 9.45 to allow for natural pH drop during procedure. If two back-to-back uses, re-pH to 9.45 before second batch of slides.

**1% AMMONIUM SULFIDE PREPARATION:** Prepare solution at time of use in fume hood. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	200uL	Concentrated Ammonium Sulfide (Not Included)	Into a chemically clean container or new/unused plasticware.
2	Add	20ml	Distilled Water	Mix thoroughly


- Carson FL Hladik C: Histotechnology A Self Instructional Text, 2009, 318 – 320
- Modifications by AMTS Research & Development Department, 1979-2019.

**STAINING PROCEDURE:** Preheating may be required. See step 3 for more information.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Immerse three separate, properly labeled serial-sectioned slides simultaneously into separate buffers as shown below.					
1a	Slide 1, Labeled: "Type 1 (4.3)"					
	Immerse	Type 1 Buffer	--	20	--	Once complete, continue to step 2.
1b	Slide 2, Labeled: "Type 1&2BC (4.6)"					
	Immerse	Type 1&2BC Buffer	--	20	--	Once complete, continue to step 2.
1c	Slide 3, Labeled: "Type 2AB (10.4)"					
	Immerse	Type 2AB Buffer	--	20	--	Once complete, continue to step 2.
2	Rinse	Running DI Water	--	1	--	
3	Immerse	ATPase Incubation Solution	37°	→	→	Slide 1 "Type 1 (4.3)" →45 Minutes. Slide 2 "Type 1&2BC (4.6)" →30 Minutes. Slide 3 "Type 2AB (10.4)" →10 Minutes. After each slide's time completes, do not rinse and continue to step 4.
4	Immerse	Calcium Chloride Solution	--	3	--	Without rinsing, continue to next step.
5	Immerse	Cobalt Chloride Solution	--	3	--	Once complete, <b>rinse in running DI water (1 minute)</b> and continue. If bath in step 2 does not change over at least once per minute, use a different bath.
6	Immerse	1% Ammonium Sulfide	--	1	--	Perform step under fume hood or exhaust-ventilated autostainer!
7	Rinse	Running Tap Water	--	1	--	
8	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
9	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
10	Coverslip	Permanent Mounting Media	--	--	--	

**Myofibrillar ATPase Procedure Flowchart**

- Carson FL Hladik C: Histotechnology A Self Instructional Text, 2009, 318 – 320
- Modifications by AMTS Research & Development Department, 1979-2019.

### AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint size or larger may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. All buffer and working solutions preparations yield approximately 20ml of solution each for the manual procedure and each must be scaled up to accommodate desired autostainer bath size. A minimum of 9 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 15 baths to run the complete procedure.

**TEST YIELD:** \*Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	225	9	250ml Glass Stain Dish	2	169	85
30ml Glass Coplin Jar	16	216	14	200ml Bath Autostainer	2	144	72
40ml Helledahl Staining Jar	12	216	18	400ml Bath Autostainer	1	153	153

### CE MARKINGS AND DESIGNATIONS:

<b>REF</b>	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
<b>LOT</b>	Batch Code		Use By	<b>EC REP</b>	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
<b>IVD</b>	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	<b>CE</b>		
	Irritant		Health Hazard			

### CONTACT INFORMATION:

#### American MasterTech Scientific

1330 Thurman St.  
Lodi, CA 95240 USA  
Tel 800 860 4073  
Fax 209 368 4136  
www.americanmastertech.com

#### StatLab

2090 Commerce Drive  
McKinney, TX 75069 USA  
Tel 800 442 3573  
Fax 972 436 1369  
www.statlab.com

1. Carson FL Hladik C: Histotechnology A Self Instructional Text, 2009, 318 – 320
2. Modifications by AMTS Research & Development Department, 1979-2019.

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**COMPOSANTS INCLUS:** Type 1 Buffer, Type 1 & 2C Buffer, Alkaline Buffer Stock, ATP Powder, Calcium Chloride Solution, Cobalt Chloride Solution, 1N Sodium Hydroxide, 1N Hydrochloric Acid. Non inclus: Ammonium Sulfide

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