Human Brain

Modified Bielschowsky's Stain Kit Procedure

100ml Kit Item #: KTBIE

Liter Kit Item#: N/A

Pint Kit Item #: N/A

Gallon Kit Item#: N/A

Control Slide(s)	Item#
Alzheimer	CSA0225P

Included Components	
Silver Nitrate Solution	Stock Developer Solution C
Stock Developer Solution A	5% Sodium Thiosulfate
Stock Developer Solution B	

Components Not Included:

Concentrated Ammonium Hydroxide (Item#: SPA1326)

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify axons, intracellular neurofibrils and senile plaques. Intracellular neurofibrils (neurofibrillary tangles), axons, and senile plaques stain gray to black, and nuclei and background yellow to gold.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 6-8µm tissue sections or frozen sections.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

AMMONIACAL SILVER SOLUTION PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	20ml	Silver Nitrate Solution	Into a chemically clean container or new/unused plasticware.
2	Add	Drops	Concentrated Ammonium Hydroxide (Not Included)	Add drop by drop while carefully agitating container until <u>no precipitate remains</u> . Solution will turn brown, then colorless.

DEVELOPER SOLUTION PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	50ml	Distilled/DI Water	Into a chemically clean container or new/unused plasticware.
2	Add	8 Drops	Stock Developer Solution "A"	Mix thoroughly.
3	Add	8 Drops	Stock Developer Solution "B"	Mix thoroughly.
4	Add	4 Drops	Stock Developer Solution "C"	Mix thoroughly.

AMMONIA WATER PREPERATION: Prepare solution at time of use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	25ml	Distilled/DI Water	Into a chemically clean container or new/unused plasticware.
2	Add	3 Drops	Concentrated Ammonium Hydroxide (Not Included)	Mix thoroughly.

STAINING PROCEDURE: Preheating Required. See steps 4 and 5 for more information.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	2	--	
4	Waterbath	Silver Nitrate Solution	40°	15	--	Immerse into preheated solution. Once complete, rinse in running DI water (2 minutes) and continue.
5	Waterbath	Ammoniacal Silver Solution	40°	10	--	Immerse into preheated solution. Once complete, drain excess solution. Without rinsing, continue to next step
6	Immerse	Developer Solution	--	--	3-30	Agitate. Remove when golden brown tissue color is achieved.
7	Immerse	Ammonia Water	--	--	30	Quickly immerse to stop developing process. Once complete, rinse in running DI water (2 minutes) and continue.
8	Immerse	5% Sodium Thiosulfate	--	2	--	Once complete, rinse in running DI water (2 minutes) and continue.
9	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
10	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minutes each change or as required if using a xylene substitute.
11	Coverslip	Permanent Mounting Media	--	--	--	

1. Mirra SS Hart MN Terry RD: Making the Diagnosis of Alzheimer's Disease, Arch Pathol Lab Med. 1993;117: 132 - 144.
2. Luna LG: Histopathologic Methods and Color Atlas of Special Stains and Tissue Artifacts; American Histolabs, Inc. 1992: 483 - 485.
3. With modifications by AMTS R&D Department, 1979-2018.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit is not available in the pint or larger size.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	12	48	4	250ml Glass Stain Dish	n/a	n/a	n/a
30ml Glass Coplin Jar	8	64	8	200ml Bath Autostainer	n/a	n/a	n/a
40ml Hellendahl Jar	6	96	16	400ml Bath Autostainer	n/a	n/a	n/a

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	StatLab Medical Products 2090 Commerce Drive McKinney, TX 75069 www.statlab.com 800-442-3573
LOT	Batch Code		Use By			
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	 GHS05	Corrosive	
 GHS07	Irritant	 GHS08	Health Hazard	 GHS09	Environmentally Damaging	

1. Mirra SS Hart MN Terry RD: Making the Diagnosis of Alzheimer's Disease, Arch Pathol Lab Med. 1993:117: 132 - 144.
2. Luna LG: Histopathologic Methods and Color Atlas of Special Stains and Tissue Artifacts; American Histolabs, Inc. 1992: 483 - 485.
3. With modifications by AMTS R&D Department, 1979-2018.

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