Human Breast

Feulgen Stain Kit Procedure

100ml Kit Item #: KTFEU

Liter Kit Item#: N/A

Pint Kit Item #: N/A

Gallon Kit Item#: N/A

Control Slide(s)
Item#

Tonsil

CST0225P

Included Components

1.0N Hydrochloric Acid

Sulfurous Rinse

Optimized Schiff's Solution™

0.2% Fast Green

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify deoxyribonucleic acid (DNA). DNA stains magenta and background light green.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-6µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

STAINING PROCEDURE: Preheating Required. See step 4 for more information.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Waterbath	1.0N Hydrochloric Acid	60°	10	--	Immerse into preheated solution. Once complete, do not rinse and continue to next step.
5	Immerse	Optimized Schiff's Solution	--	15	--	Without rinsing, continue to next step.
6	Rinse	Sulfurous Rinse, 2 changes	--	3	--	3 minutes each change.
7	Rinse	Running Tap Water	--	5	--	
8	Immerse	0.2% Fast Green	--	--	10	Without rinsing, continue to next step.
9	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	Fresh changes. 1 minute each change.
10	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
11	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB; Theory and practice of Histotechnology; 1980, 150.
2. With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit is not available in the pint or larger size.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	225	9	250ml Glass Stain Dish	n/a	n/a	n/a
30ml Glass Coplin Jar	16	216	14	200ml Bath Autostainer	n/a	n/a	n/a
40ml Hellendahl Jar	12	216	18	400ml Bath Autostainer	n/a	n/a	n/a

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		
 GHS02	Flammable	 GHS08	Health Hazard			

CONTACT INFORMATION:

American MasterTech Scientific
 1330 Thurman St.
 Lodi, CA 95240 USA
 Tel 800 860 4073
 Fax 209 368 4136
 www.americanmastertech.com

StatLab
 2090 Commerce Drive
 McKinney, TX 75069 USA
 Tel 800 442 3573
 Fax 972 436 1369
 www.statlab.com

1. Sheehan DC Hrapchak BB; Theory and practice of Histotechnology; 1980, 150.
2. With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: 1.0N Hydrochloric Acid, Optimized Schiff's Solution™, Sulfurous Rinse, 0.2% Fast Green

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