Human Colon

G.M.S. Stain Kit Procedure

100ml Kit Item #: KTGMS

Liter Kit Item#: KTGMSLT

Pint Kit Item #: KTGMSPT

Gallon Kit Item#: N/A

Control Slide(s)
Item#

Aspergillus	CSA0625P
Candida Albicans	CSC0525P
Pneumocystis	CSP0125P

Included Components

Methenamine/Silver Compound	0.3% Gold Chloride
Sodium Borate Concentrate	2% Sodium Thiosulfate
5% Chromic Acid	Light Green Counterstain
1% Sodium Bisulfite	Slide Jar

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify fungi and *Pneumocystis jirovecii*. Fungi and *Pneumocystis jirovecii* stain gray to black and background green.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

METHENAMINE/SILVER SOLUTION PREPARATION: Prepare solution at time of use. Solution expires after one use. Methenamine/Silver Compound may soften during transit; this will not affect product performance. If desired, Methenamine/Silver Compound can be re-hardened by placing vials in a freezer. Methenamine/Silver Solution reacts with metals. Do not use metal forceps, racks, or tools with the solution.

#	Action	Amount	Chemical/Reagent	Details
1	Add	2ml	Distilled Water	Into a vial of Methenamine/Silver Compound. Break up chunks with an applicator stick. Pour contents of the vial into a chemically cleaned container or new/unused plasticware.
2	Add	46ml	Distilled Water	Mix thoroughly.
3	Add	2ml	Sodium Borate Concentrate	Mix for 2 minutes.

STAINING PROCEDURE: Preheating Required. See step 6 for more information. Microwave procedure located on opposite page.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Immerse	5% Chromic Acid	--	60	--	Once complete, rinse in running DI water (1 minute) and continue.
5	Immerse	1% Sodium Bisulfite	--	1	--	Once complete, rinse in running DI water (1 minute) and continue.
6	Waterbath	Methenamine/Silver Solution	62°	15-25	--	Immerse into preheated solution. Once complete, rinse in running DI water (1 minute) and continue.
7	Immerse	0.3% Gold Chloride	--	1	--	Once complete, rinse in running DI water (1 minute) and continue.
8	Immerse	2% Sodium Thiosulfate	--	1	--	Once complete, rinse in running DI water (1 minute) and continue.
9	Immerse	Light Green Counterstain	--	--	30-60	
10	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	Fresh changes. 1 minute each change.
11	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
12	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB; Theory and Practice of Histotechnology; 1980, 245 - 246.
2. Sale GE: Rapid Methenamine Silver Stain; Arch Pathol Lab Med: 102; 1978, 351 - 352.
3. With modifications by AMTS R&D Department, 1979-2019.

MICROWAVE STAINING PROCEDURE: Stain no more than 3 slides at a time when using the supplied slide jars. To allow sufficient volume for thermal expansion, do not fill slide jars past the bottom cap thread.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Microwave	5% Chromic Acid	70-80°	--	--	Immerse slide into the supplied slide jar containing 5% Chromic Acid, loosely cap, and heat (do not boil). Incubate for 1 minute. Once complete, rinse in running DI water (1 minute) and continue.
5	Immerse	1% Sodium Bisulfite	--	1	--	Once complete, rinse in running DI water (1 minute) and continue.
6	Microwave	Methenamine/Silver Solution	65-75°	--	--	Immerse slide into the supplied slide jar containing Methenamine/Silver Solution, loosely cap, and heat (do not boil). Incubate for 5 minutes (maintain temperature) until tissue is medium brown. Once complete, rinse in running DI water (1 minute) and continue.
7	Immerse	0.3% Gold Chloride	--	1	--	Once complete, rinse in running DI water (1 minute) and continue.
8	Immerse	2% Sodium Thiosulfate	--	1	--	Once complete, rinse in running DI water (1 minute) and continue.
9	Immerse	Light Green Counterstain	--	--	30-60	
10	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	Fresh changes. 1 minute each change.
11	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
12	Coverslip	Permanent Mounting Media	--	--	--	

AUTOSTAINER CONFIGURATION AND NOTES:

This kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the procedure grid on the reverse side of this page. Standard Methenamine/Silver Solution preparation procedure yields 50ml of solution and must be scaled up to accommodate desired autostainer bath size. A minimum of 7 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 18 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	50	200	4	250ml Glass Stain Dish	4	120	30
30ml Glass Coplin Jar	33	264	8	200ml Bath Autostainer	5	150	30
40ml Hellendahl Jar	25	400	16	400ml Bath Autostainer	2	100	50

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		
	Oxidizing		Irritant		Health Hazard	
GHS03		GHS07		GHS08		

CONTACT INFORMATION:

American MasterTech Scientific
 1330 Thurman St. • Lodi, CA 95240 USA
 Tel 800 860 4073 • Fax 209 368 4136
www.americanmastertech.com

StatLab
 2090 Commerce Drive • McKinney, TX 75069 USA
 Tel 800 442 3573 • Fax 972 436 1369
www.statlab.com

1. Sheehan DC Hrapchak BB; Theory and Practice of Histotechnology; 1980, 245 - 246.
2. Sale GE: Rapid Methenamine Silver Stain; Arch Pathol Lab Med: 102; 1978, 351 - 352.
3. With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: Methenamine/Silver Compound, Sodium Borate Concentrate, 5% Chromic Acid, 1% Sodium Bisulfate, 0.3% Gold Chloride, 2% Sodium Thiosulfate, Light Green Counterstain, Jarres de diapositives.

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