Human Lung

Modified Gridley's Stain Kit Procedure

100ml Kit Item #: KTGRI

Liter Kit Item#: N/A

Pint Kit Item #: N/A

Gallon Kit Item#: N/A

Control Slide(s)
Item#

Amoeba	CSA0325P
Aspergillus	CSA0625P
Candida Albicans	CSC0525P

Included Components

1% Periodic Acid Solution	Aldehyde Fuchsin Solution
Optimized Schiff's Solution	Metanil Yellow Solution
Sulfurous Rinse	

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify hyphae, conidia, elastin and mucin. Hyphae are stained magenta to purple, conidia dark rose to purple, elastic tissue and mucin purple, and background yellow.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 5-6µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Immerse	1% Periodic Acid Solution	--	30	--	Once complete, rinse in running DI water (3 minutes) and continue.
5	Immerse	Optimized Schiff's Solution	--	15	--	
6	Rinse	Sulfurous Rinse, 2 changes	--	2	--	2 minutes each change. Once complete, rinse in running DI water (5 minutes) and continue.
7	Immerse	Aldehyde Fuchsin Solution	--	15-20	--	
8	Rinse	95% Reagent Alcohol	--	--	10-30	To remove excess stain. Once complete, rinse in running DI water until no stain runs off slide and continue.
9	Immerse	Mentanil Yellow Solution	--	--	30-60	Once complete, rinse in running DI water (1 minute) and continue.
10	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
11	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
12	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB; Theory and Practice of Histotechnology; 1980, 245 - 246.
2. With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit is not available in the pint or larger size.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	5	100	20	250ml Glass Stain Dish	n/a	n/a	n/a
30ml Glass Coplin Jar	3	90	30	200ml Bath Autostainer	n/a	n/a	n/a
40ml Hellendahl Jar	2	80	40	400ml Bath Autostainer	n/a	n/a	n/a

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		
 GHS07	Irritant	 GHS08	Health Hazard			

CONTACT INFORMATION:

American MasterTech Scientific
 1330 Thurman St.
 Lodi, CA 95240 USA
 Tel 800 860 4073
 Fax 209 368 4136
 www.americanmastertech.com

StatLab
 2090 Commerce Drive
 McKinney, TX 75069 USAW
 Tel 800 442 3573
 Fax 972 436 1369
 www.statlab.com

1. Sheehan DC Hrapchak BB; Theory and Practice of Histotechnology; 1980, 245 - 246.
2. With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: 0.1% Periodic Acid Solution, Optimized Schiff's Solution, Sulfurous Rinse, Aldehyde Fuchsin Solution, Metanil Yellow Solution.

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COMPONENTES INCLUIDOS: 0.1% Periodic Acid Solution, Optimized Schiff's Solution, Sulfurous Rinse, Aldehyde Fuchsin Solution, Metanil Yellow Solution.

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