Human Skeletal Muscle

Masson's Trichrome Stain Kit Procedure

100ml Kit Item #: KTMTR
Pint Kit Item #: KTMTRPT

Liter Kit Item#: KTMTRLT
Gallon Kit Item#: N/A

Control Slide(s)	Item#	Included Components	
Heart	CSH0725P	Weigert's Hematoxylin A	Phosphomolybdic /
Skeletal Muscle	CSS0725P	Weigert's Hematoxylin B	Phosphotungstic Acid
Trichrome	CST0125P	Bouin's Fluid	Aniline Blue Stain
Uterus	CSU0325P	Biebrich Scarlet - Acid Fuchsin	1% Acetic Acid

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify collagen and muscle. Cytoplasm, keratin, muscle, and intercellular fiber are stained red, collagen and mucus blue, and nuclei black.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-6µm tissue section

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

WORKING WEIGERT'S HEMATOXYLIN PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	25ml	Weigert's Hematoxylin A	Into a chemically clean container or new/unused plasticware.
2	Add	25ml	Weigert's Hematoxylin B	Mix thoroughly.

STAINING PROCEDURE: Preheating required. See step 4 for more information. Microwave procedure located on opposite page.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	Waterbath	Bouin's Fluid	56-62°	60	--	Immerse into preheated solution. (Optional: Immerse overnight at room temperature). Once complete, rinse in running tap water (3 minutes) until tissue is colorless and continue.
5	Immerse	Working Weigert's Hematoxylin	--	5	--	Once complete, rinse in running tap water (2 minutes) and continue.
6	Immerse	Biebrich Scarlet-Acid Fuchsin	--	15	--	If staining Central Nervous System (C.N.S. Sections), reduce staining time to 1-2 minutes. Once complete, rinse in running tap water (1 minute) and continue.
7	Immerse	Phosphomolybdic / Phosphotungstic Acid	--	10-15	--	Without rinsing continue to next step.
8	Immerse	Aniline Blue Stain	--	5-10	--	If staining Central Nervous System (C.N.S. Sections), increase staining time to 15-20 minutes. Once complete, rinse in running tap water (1 minute) and continue.
9	Immerse	1% Acetic Acid	--	3-5	--	
10	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	Fresh changes. 1 minute each change.
11	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
12	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 190.
2. A.F.I.P. Laboratory Methods in Histotechnology: 1992, 132 - 133.
3. With modifications by AMTS R&D Department, 1979-2019.

MICROWAVE STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	Microwave	Bouin's Fluid	80-85°	--	--	Heat Bouin's Fluid in a loosely capped slide jar (without slides). Do not boil. Remove from microwave, immerse slide, tighten cap, agitate, and incubate for 4 minutes. Once complete, rinse in running tap water (3 minutes) until tissue is colorless and continue.
5	Microwave	Working Weigert's Hematoxylin	80-90°	--	--	Immerse slide into loosely capped jar of Weigert's Hematoxylin and heat (do not boil). Remove from microwave, cap jar tightly, and incubate for 30 seconds. Once complete, rinse in running tap water (1 minute) and continue.
6	Microwave	Biebrich Scarlet-Acid Fuchsin	80-90°	--	--	Immerse slide into loosely capped jar of Biebrich Scarlet-Acid Fuchsin and heat (do not boil). Remove from microwave, cap jar tightly, and incubate for 2 minutes. Once complete, rinse in running tap water (1 minute) and continue.
7	Microwave	Phosphomolybdic / Phosphotungstic Acid	80-90°	--	--	Immerse slide into loosely capped jar of Phosphomolybdic / Phosphotungstic Acid and heat (do not boil). Remove from microwave, cap jar tightly, and incubate for 1 minute. Without rinsing continue to step 8.
8	Microwave	Aniline Blue Stain	80-90°	--	--	Immerse slide into loosely capped jar of Aniline Blue Stain and heat (do not boil). Remove from microwave, cap jar tightly, and incubate for 45-90 seconds. Once complete, rinse in running tap water (1 minute) and continue.
9	Immerse	1% Acetic Acid (Room Temperature)	--	3-5	--	
10	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	Fresh changes. 1 minute each change.
11	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
12	Coverslip	Permanent Mounting Media	--	--	--	

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. A minimum of 7 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 18 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	500	20	250ml Glass Stain Dish	2	375	188
30ml Glass Coplin Jar	16	480	30	200ml Bath Autostainer	2	320	160
40ml Hellendahl Jar	12	480	40	400ml Bath Autostainer	1	340	340

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use		Flammable	
	Corrosive		Irritant		Health Hazard	

CONTACT INFORMATION:

American MasterTech Scientific
 1330 Thurman St. • Lodi, CA 95240 USA
 Tel 800 860 4073 • Fax 209 368 4136 • www.americanmastertech.com

StatLab
 2090 Commerce Drive • McKinney, TX 75069 USA
 Tel 800 442 3573 • Fax 972 436 1369 • www.statlab.com

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 190.
2. A.F.I.P. Laboratory Methods in Histotechnology: 1992, 132 - 133.
3. With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: Weigert's Hematoxylin A, Weigert's Hematoxylin B, Bouin's Fluid, Biebrich Scarlet – Acid Fuchsin, Phosphomolybdic/Phosphotungstic Acid, Aniline Blue Stain, 1% Acetic Acid

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
1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 190.
2. A.F.I.P. Laboratory Methods in Histotechnology; 1992, 132 - 133.
3. With modifications by AMTS R&D Department, 1979-2019.

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COMPONENTES INCLUIDOS: Weigert's Hematoxylin A, Weigert's Hematoxylin B, Bouin's Fluid, Biebrich Scarlet – Acid Fuchsin, Phosphomolybdic/Phosphotungstic Acid, Aniline Blue Stain, 1% Acetic Acid

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