Human Uterus

P.T.A.H. Stain Kit Procedure, Microwave

100ml Kit Item #: KTPTA

Liter Kit Item#: N/A

Pint Kit Item #: KTPTAPT

Gallon Kit Item#: N/A

Control Slide(s)	Item#
Heart	CSH0725P
Skeletal Muscle	CSS0725P
Trichrome	CST0125P
Uterus	CSU0325P

Included Components
10% Zinc Chloride Solution
5% Ferric Ammonium Sulfate Solution
P.T.A.H. Stain


PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify collagen, fibrin and striated muscle. Fibrin and muscle stain blue, collagen brownish-red to red, and nuclei blue.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

MICROWAVE STAINING PROCEDURE: Standard staining procedure located on opposite page.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Microwave	10% Zinc Chloride Solution	80-90°	--	--	Heat 10% Zinc Chloride Solution in a loosely capped slide jar (without slide). Do not boil. Remove from microwave, tighten cap, and agitate to equalize temperature. Immerse slide and incubate for 15 minutes while agitating periodically. Once complete, rinse in running DI water (1 minute) and continue.
5	Microwave	5% Ferric Ammonium Sulfate Solution	80-90°	--	--	Immerse slide into a loosely capped jar of 5% Ferric Ammonium Sulfate Solution and heat (do not boil). Incubate for 2 minutes. Once complete, rinse in running DI water (1 minute) and continue.
6	Microwave	P.T.A.H. Stain	80-90°	--	--	Heat P.T.A.H. Stain in a loosely capped slide jar (without slide). Do not boil. Remove from microwave, tighten cap, and agitate to equalize temperature. Immerse slide and incubate for 15 minutes while agitating periodically. Re-heat solution (with slide) to 80-90°C (do not boil). Remove from microwave, agitate to equalize temperature and incubate for 15 more minutes. Once complete, rinse in running DI water (1 minute) and continue.
7	Differentiate	95% Reagent Alcohol	--	--	30-60	Check using microscope for correct differentiation.
8	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
9	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
10	Coverslip	Permanent Mounting Media	--	--	--	

- Shapiro SH Sohn LC: Rapid Microwave Phosphotungstic Acid-Hematoxylin Stain for Paraffin and Glycol Methacrylate Sections; The Journal of Histotechnology; Vol 17, No 2, June 1994; 125-126.
- A.F.I.P. Laboratory Methods in Histotechnology: 1992; 98 - 99.
- With modifications by AMTS R&D Department, 1979-2019.

STANDARD STAINING PROCEDURE: Preheating Required. See step 6 for more information.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohols, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Immerse	10% Zinc Chloride Solution	--	180	--	Once complete, rinse in running DI water (1 minute) and continue.
5	Immerse	5% Ferric Ammonium Sulfate Solution	--	10	--	Without rinsing, continue to next step.
6	Waterbath	P.T.A.H. Stain	65°	120	--	Immerse into preheated solution. For room temperature staining without a waterbath, increase time to 12-24 hours while periodically agitating.
7	Differentiate	95% Reagents Alcohol	--	--	30-60	Check using microscope for correct differentiation.
8	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
9	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
10	Coverslip	Permanent Mounting Media	--	--	--	

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on this page. A minimum of 5 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 15 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	225	9	250ml Glass Stain Dish	2	169	85
30ml Glass Coplin Jar	16	216	14	200ml Bath Autostainer	2	144	72
40ml Hellendahl Jar	12	216	18	400ml Bath Autostainer	1	153	153

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		
	Corrosive		Environmentally Damaging			
GHS05		GHS09				

CONTACT INFORMATION:

American MasterTech Scientific
 1330 Thurman St.
 Lodi, CA 95240 USA
 Tel 800 860 4073
 Fax 209 368 4136
 www.americanmastertech.com

StatLab
 2090 Commerce Drive
 McKinney, TX 75069 USA
 Tel 800 442 3573
 Fax 972 436 1369
 www.statlab.com

- Shapiro SH Sohn LC: Rapid Microwave Phosphotungstic Acid-Hematoxylin Stain for Paraffin and Glycol Methacrylate Sections; The Journal of Histotechnology; Vol 17, No 2, June 1994; 125-126.
- A.F.I.P. Laboratory Methods in Histotechnology: 1992; 98 - 99.
- With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: 10% Zinc Chloride Solution, 5% Ferric Ammonium Sulfate Solution, P.T.A.H. Stain

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