Human Colon Tissue

Russell-Movat Pentachrome Stain Kit Procedure

100ml Kit Item #: KTRMP

Liter Kit Item#: N/A

Pint Kit Item #: KTRMPPT

Gallon Kit Item#: N/A

Control Slide(s)	Item#	Included Components	
Aorta	CSA0825P	10% Alcoholic Hematoxylin	3% Acetic Acid
Skeletal Muscle	CSS0725P	2% Ferric Chloride	1% Alcian Blue Solution
Normal Colon	CSC1025P	10% Ferric Chloride	Crocein Scarlet-Acid Fuchsin
		Universal Iodine Solution™	5% Phosphotungstic Acid
		5% Sodium Thiosulfate	Alcoholic Saffron Solution
		1% Acetic Acid	

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify collagen, elastic fibers, fibrinoid, mucin, muscle, and nuclei. Collagen is stained yellow, elastic fibers and nuclei black, fibrinoid intense red, mucins blue to green, and muscle red.


SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

VERHOEFF'S ELASTIC STAIN PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	10ml	10% Alcoholic Hematoxylin	Into a chemically clean container or new/unused plasticware.
2	Add	10ml	Reagent Alcohol, <i>Not Included</i>	Mix thoroughly.
3	Add	10ml	10% Ferric Chloride	Mix thoroughly.
4	Add	10ml	Universal Iodine Solution™	Mix thoroughly.
5	Filter	--	Paper Towel	Filter into a chemically cleaned container or new/unused plasticware.

STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Immerse	Verhoeff's Elastic Stain	--	15	--	Once complete, rinse in running DI water (2 minutes) and continue.
5	Immerse	2% Ferric Chloride	--	--	60	Until elastic fibers are sharply defined (use microscope). Once complete, rinse in running DI water (1 minute) and continue.
6	Immerse	5% Sodium Thiosulfate	--	1	--	Once complete, rinse in running DI water (2 minutes) and continue.
7	Immerse	3% Acetic Acid	--	3	--	Without rinsing continue to next step.
8	Immerse	1% Alcian Blue Solution	--	15	--	Until mucins are blue. Once complete, rinse in running DI water (1 minute) and continue.
9	Immerse	Crocein Scarlet-Acid Fuchsin	--	2	--	Once complete, rinse in running DI water (1 minute) and continue.
10	Immerse	1% Acetic Acid	--	--	1-2	5 dips for 1-2 seconds each.
11	Immerse	5% Phosphotungstic Acid, 2 changes	--	2	--	2 minutes each change. Check section under microscope and stop differentiation when connective tissue is clear, before the elastic fibers are de-stained. Once complete, immerse in 1% Acetic Acid (5 dips for 1-2 seconds each) and continue.
12	Immerse	Absolute Alcohol, 2 changes	--	1	--	1 minute each change.
13	Immerse	Alcoholic Saffron Solution	--	15	--	
14	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	Fresh changes. 1 minute each change.
15	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
16	Coverslip	Permanent Mounting Media	--	--	--	

1. Movat HZ: Demonstration of all connective tissue elements in a single section; Arch Path 60: 1955, 289.
2. Russell K: Pentachrome stain modification; Arch Path 94: 1972, 187.
3. With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. Standard Verhoeff's Elastic Stain preparation procedure yields approximately 40ml of solution and must be scaled up to accommodate desired autostainer bath size. A minimum of 13 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 24 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	225	9	250ml Glass Stain Dish	2	169	85
30ml Glass Coplin Jar	16	216	14	200ml Bath Autostainer	2	144	72
40ml Helledahl Staining Jar	12	216	18	400ml Bath Autostainer	1	153	153

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use		Flammable	CE
	Corrosive		Irritant		Health Hazard	

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. Movat HZ: Demonstration of all connective tissue elements in a single section; Arch Path 60: 1955, 289.
2. Russell K: Pentachrome stain modification; Arch Path 94: 1972, 187.
3. With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: 10% Alcoholic Hematoxylin, Universal Iodine Solution™, 5% Sodium Thiosulfate, 1% Alcian Blue Solution, 1% Acetic Acid, Alcoholic Saffron Solution, 10% Ferric Chloride, 3% Acetic Acid, 2% Ferric Chloride, Crocein Scarlet-Acid Fuchsin, Phosphotungstic Acid

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