American MasterTech
scientific laboratory supplies

AquaSlip™

Aqueous Permanent Mounting Medium

(Rev. 3/8/18)

AquaSlip™, The product and its purpose:

AEC and Fast Red are commonly used chromogens for peroxidase and alkaline-phosphatase based IHC immunostain systems. However, slides stained with these chromogens cannot be coverslipped using organic solvent based mounting medium. AquaSlip™ by American MasterTech is an aqueous mounting medium with an optimal refractive index, which, when applied to the stained tissue sections, permanently protects the section without the need for a coverslip. If desired, a coverslip can be applied following complete drying of AquaSlip™ by using an organic solvent based permanent mounting medium. Whether coverslipped or not, AquaSlip™ provides crisp, clear, optimum appearance of stained tissue. AquaSlip™ is compatible with DAB, AEC, BCIP/INT, BCIP/NBT, and Fast Red. AquaSlip™ is designed for permanent mounting of tissue specimens stained using peroxidase or alkaline-phosphatase based detection systems. AquaSlip™ is intended for In Vitro Diagnostic use. Active ingredients include; Homopolymer.

AquaSlip™, Instructions for use:

Mounting Procedure:

Place AquaSlip™ bottle upside down in a container before using, which will clear AquaSlip™ of tiny bubbles. Blot excess water from the glass without letting the tissue sections dry; insure that tissue section remains wet prior to applying AquaSlip™. Apply 3 drops of AquaSlip™ on tissue section, and then gently rotate the slide using wrist movement to form a thin layer of AquaSlip™; insure that the tissue is completely covered. Do not apply coverslip. Place slides in a 70° C preheated oven. NOTE: Prior to preheating oven, insure that the oven is level. Dry slides for 20 minutes, then check AquaSlip™ for complete drying. Return slides to oven if more drying time is necessary. Remove slides from oven. Allow them to cool to room temperature (RT). Evaluate sections using a microscope.

Post Mounting:

Once AquaSlip™ has thoroughly dried, apply an organic based mounting media such as MasterTech ClearMount™, CoverSafe™, or MicroMount™ and then apply coverslip. Because the refractive index of AquaSlip™ is very close to organic solvent based mounting mediums, no loss of clarity occurs when post mounting.

Removal of AquaSlip™:

AquaSlip™ is an aqueous mounting medium and can be easily removed by placing the slides in a beaker or coplin jar filled with deionized water for 1 to 2 hours, then move slides to a fresh container of deionized water and agitate slides gently one at a time to complete removal of AquaSlip™.

AquaSlip™, Storage and stability:

Store AquaSlip™ at 15°C to 30°C. Do not refrigerate!

AquaSlip™, Precautions:

See Safety Data Sheet (SDS) for warnings and precautions.

AquaSlip™, Disposal:

Dispose of AquaSlip™ in accordance with all local, state, federal, environmental and/or EU regulations.

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CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use			

Copyright © 2015-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.