Normal Human Colon

## Alcian Blue pH 2.5 Stain Procedure

**100ml Item#:** STABL2.5100

**Liter Item#:** STABL2.5LT

**Pint Item#:** STABL2.5PT

**Gallon Item#:** STABL2.5GAL

| Control Slide(s) | Item#    | Required Components (Not Included) |
|------------------|----------|------------------------------------|
| Colon            | CSC1025P | Glacial Acetic Acid                |
| Small Intestine  | CSS0325P | Nuclear Fast Red Stain             |

**PRINCIPLE AND RESULTS:** This stain is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify acid mucins, carboxylated and sulfated mucosubstances. Strongly sulfated mucosubstances are not demonstrated. Acid mucins, carboxylated and sulfated mucosubstances are stained blue.

**SPECIMEN CRITERIA:** Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

**STORAGE AND USAGE NOTES:** Store/Use according to the temperature and expiration on the label.

**PRECAUTIONS:** For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

**3% ACETIC ACID PREPARATION:** Prepare solution at time of use. Solution expires after one use.

| # | Action | Amount | Chemical/Reagent    | Details                                                        |
|---|--------|--------|---------------------|----------------------------------------------------------------|
| 1 | Add    | 97ml   | Distilled Water     | Into a chemically cleaned container or new/unused plasticware. |
| 2 | Add    | 3ml    | Glacial Acetic Acid | Mix thoroughly.                                                |

### STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

| # | Action        | With                            | Heat °C | Time |      | Details                                                                   |
|---|---------------|---------------------------------|---------|------|------|---------------------------------------------------------------------------|
|   |               |                                 |         | Mins | Secs |                                                                           |
| 1 | Deparaffinize | Xylene or Substitute, 2 changes | --      | 5    | --   | 5 minutes each change or as required if using a xylene substitute.        |
| 2 | Rinse         | Absolute Alcohol, 3 changes     | --      | 1    | --   | 1 minute each change or as required if using graded alcohols.             |
| 3 | Rinse         | Running Tap Water               | --      | 1    | --   |                                                                           |
| 4 | Immerse       | 3% Acetic Acid                  | --      | 3    | --   | Without rinsing, continue to next step.                                   |
| 5 | Immerse       | Alcian Blue Stain pH 2.5        | --      | 30   | --   | Once complete, <b>rinse in running tap water (1 minute)</b> and continue. |
| 6 | Immerse       | Nuclear Fast Red Stain          | --      | 5    | --   | Once complete, <b>rinse in running tap water (1 minute)</b> and continue. |
| 7 | Dehydrate     | Absolute Alcohol, 3 changes     | --      | 1    | --   | 1 minute each change.                                                     |
| 8 | Clear         | Xylene or Substitute, 3 changes | --      | 1    | --   | 1 minute each change or as required if using a xylene substitute.         |
| 9 | Coverslip     | Permanent Mounting Media        | --      | --   | --   |                                                                           |

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology, 1980, 172 – 173
2. With modifications by AMTS R&D Department, 1979-2019.

**CE MARKINGS AND DESIGNATIONS:**

|            |                                    |                                                                                   |                                   |                                                                                            |                |                                                                                                                  |
|------------|------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------|--------------------------------------------------------------------------------------------|----------------|------------------------------------------------------------------------------------------------------------------|
| <b>REF</b> | Catalogue Number                   |  | Temperature Limitation            |           | Manufacturer   | American MasterTech Scientific<br>1330 Thurman St.<br>Lodi, CA 95240 USA<br>Tel 800 860 4073<br>Fax 209 368 4136 |
| <b>LOT</b> | Batch Code                         |  | Use By                            | <b>EC REP</b>                                                                              | Representative | Emergo Europe<br>Prinsessegracht 20<br>2514 AP The Hague<br>The Netherlands                                      |
| <b>IVD</b> | In Vitro Diagnostic Medical Device |  | Consult Instructions Prior to Use | <br>GHS08 | Health Hazard  | <b>CE</b>                                                                                                        |

**CONTACT INFORMATION:**

**American MasterTech Scientific**  
 1330 Thurman St.  
 Lodi, CA 95240 USA  
 Tel 800 860 4073  
 Fax 209 368 4136  
[www.americanmastertech.com](http://www.americanmastertech.com)

**StatLab**  
 2090 Commerce Drive  
 McKinney, TX 75069 USA  
 Tel 800 442 3573  
 Fax 972 436 1369  
[www.statlab.com](http://www.statlab.com)

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology, 1980, 172 – 173
2. With modifications by AMTS R&D Department, 1979-2019.

## MULTILANGUAGE PROCEDURE

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**COMPOSANTS NON INCLUS:** Acide acétique glacial, Nuclear Fast Red Stain

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