Human Breast

Aldehyde Fuchsin Stain Procedure

100ml Item #: STAFU100

Liter Item#: N/A

Pint Item #: N/A

Gallon Item#: N/A

Control Slide(s)
Item#

Aorta

CSA0825P

Mast Cell

CSM0125P

Required Components (Not Included)

70% Reagent Alcohol

Light Green Counterstain

PRINCIPLE AND RESULTS: This stain is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify elastic fibers, beta cells, mast cell granules and pancreatic islet cells. Elastic Fibers, Mast Cell granules, Beta Cells, and Pancreatic Islet Cells are stained violet.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 5-6µm tissue section.

STORAGE AND USAGE NOTES: Store/Use component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Immerse	Aldehyde Fuchsin Stain	--	15	--	Beta cells and Pancreatic Islet cells require 30 minutes to 2 hours.
4	Rinse	 Running DI Water	--	--	5-10	Check for desired intensity (use microscope). Repeat step if necessary.
5	Rinse	70% Reagent Alcohol	--	--	5-10	Rinse quickly to remove excess stain. Once complete, rinse in running DI Water (10 seconds) and check for desired intensity (use microscope). If deeper intensity is desired, rinse in 70% Reagent Alcohol, return to Aldehyde Fuchsin Stain, and repeat steps 4-5.
6	Immerse	Light Green Counterstain	--	--	30-60	
7	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
8	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
9	Coverslip	Permanent Mounting Media	--	--	--	

- Gomori G: Aldehyde Fuchsin - A New Stain for Elastic Fibers; Am J Clin Pathol; 1950, 20: 665 - 666.
- With modifications by AMTS R&D Department, 1979-2019.

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	 GHS08	Health Hazard	

CONTACT INFORMATION:

American MasterTech Scientific
 1330 Thurman St.
 Lodi, CA 95240 USA
 Tel 800 860 4073
 Fax 209 368 4136
www.americanmastertech.com

StatLab
 2090 Commerce Drive
 McKinney, TX 75069 USA
 Tel 800 442 3573
 Fax 972 436 1369
www.statlab.com

1. Gomori G: Aldehyde Fuchsin - A New Stain for Elastic Fibers; Am J Clin Pathol; 1950, 20: 665 - 666.
2. With modifications by AMTS R&D Department, 1979-2019.

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