Chondroma

0.2% Fast Green Stain Procedure

0:05

100ml Item#: STFGS100

Liter Item#: N/A

Pint Item#: STFGSPT

Gallon Item#: N/A

Control Slide(s)
Item#

Chondroma

CSC1125P

Required Components (Not Included)

Weigert's Hematoxylin A

1% Safranin O

Weigert's Hematoxylin B

1% Acetic Acid

PRINCIPLE AND RESULTS: This stain is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to stain cytoplasm and background tissue green in the MasterTech Safranin O Stain Kit (Item#: KTSFO).

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

WEIGERT'S HEMATOXYLIN PREPARATION: Prepare solution at time of use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	1/1 Ratio	Weigert's Hematoxylin A	Into a chemically clean container or new/unused plasticware.
2	Add	1/1 Ratio	Weigert's Hematoxylin B	Mix thoroughly.

STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	Immerse	Weigert's Hematoxylin	--	5	--	Once complete, rinse in running tap water (1 minute) and continue.
5	Immerse	0.2% Fast Green	--	5	--	Without rinsing, continue to next step.
6	Immerse	1% Acetic Acid	--	--	10-15	Agitate. Without rinsing, continue to next step.
7	Immerse	1% Safranin O	--	5	--	
8	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
9	Clear	Xylene/Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
10	Coverslip	Permanent Mounting Media	--	--	--	

- Carson FL Hladik C: Histotechnology, A Self-Instructional Text; 3rd Edition, 2009, 236.
- With modifications by AMTS R&D Department, 1979-2019.

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use			
 GHS02	Flammable	 GHS08	Health Hazard			

CONTACT INFORMATION:

American MasterTech Scientific
 1330 Thurman St.
 Lodi, CA 95240 USA
 Tel 800 860 4073
 Fax 209 368 4136
www.americanmastertech.com

StatLab
 2090 Commerce Drive
 McKinney, TX 75069 USA
 Tel 800 442 3573
 Fax 972 436 1369
www.statlab.com

1. Carson FL Hladik C: Histotechnology, A Self-Instructional Text; 3rd Edition, 2009, 236.
2. With modifications by AMTS R&D Department, 1979-2019.

MULTILINGUE PROCEDURE

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COMPOSANTS NON INCLUS: Weigert's Hematoxylin A, Weigert's Hematoxylin B, 1% Safranin o, 1% Acetic Acid.

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