Human Bone Marrow

Giemsa Stain (Stock) Procedure

100ml Item #: STGIE100

Liter Item#: N/A

Pint Item #: STGIEPT

Gallon Item#: N/A

Control Slide(s)
Item#

Bone Marrow CSB0325P

Helicobacter Pylori CSH0125P

Mast Cell CSM0125P

Toxoplasma Gondii CST0325P

Required Components (Not Included)

Phosphate Buffer pH 6.8

May-Grunwald Stock

 0.25% Acetic Acid (*Mast Cells only*)

PRINCIPLE AND RESULTS: This stain is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify bone marrow cells, *Helicobacter pylori*, mast cells, rickettsias, and *Toxoplasma gondii*. When used in the May-Grunwald – Giemsa procedure, Nuclei is stained blue to violet, rickettsias intense reddish-purple, cytoplasm light blue, collagen and muscle pale pink, erythrocytes gray, yellow or pink, *Helicobacter pylori* blue, mast cells dark blue with red granules, and *Toxoplasma gondii* blue.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 5µm tissue section, or air dried smear.

STORAGE AND USAGE NOTES: Store/Use component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

WORKING MAY-GRUNWALD STAIN PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	25ml	May-Grunwald Stock	Into a coplin jar or container.
2	Add	25ml	Phosphate Buffer pH 6.8	Mix thoroughly.

WORKING GIEMSA STAIN PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	50ml	Phosphate Buffer pH 6.8	Into a coplin jar or container.
2	Add	2.5ml	Giemsa Stock	Mix thoroughly.

STAINING PROCEDURE: If staining a smear, begin procedure at step #3.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Immerse	Working May-Grunwald Stain	--	6	--	Or flood on slide rack. Agitate every few minutes for optimum staining.
5	Rinse	Phosphate Buffer pH 6.8	--	--	5-10	Or flood on slide rack until no stain runs off slide.
6	Immerse	Working Giemsa Stain	--	13	--	Or flood on slide rack. Agitate every few minutes for optimum staining.
7	Rinse	Phosphate Buffer pH 6.8	--	--	5-10	Or flood on slide rack until no stain runs off slide.
8	To stain for mast cells, follow procedure (8a) below. For all other staining, skip to step 9.					

8a. Mast Cell Procedure

Immerse	0.25% Acetic Acid	--	--	1	3-5 dips for 1 second each until desired background intensity is achieved. Use a quick Phosphate Buffer pH 6.8 rinse and microscope to check for correct differentiation and continue at step 9.
---------	-------------------	----	----	---	--

9	Immerse	Phosphate Buffer pH 6.8	--	3	--	Or flood on slide rack.
10	Air Dry	On Benchtop, or in Empty Bath	--	5-10	--	Until section is completely dry.
11	Clear	Xylene or Substitute	--	--	1	5-10 dips for 1 second each.
12	Coverslip	Permanent Mounting Media	--	--	--	

1. A.F.I.P. Manual: pg 111. Strumia: Laboratory and Clinical Medicine; 930.
2. With modifications by AMTS R&D Department, 1979-2019.

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		
	Flammable		Toxic		Health Hazard	
GHS02		GHS06		GHS08		

CONTACT INFORMATION:

American MasterTech Scientific
 1330 Thurman St.
 Lodi, CA 95240 USA
 Tel 800 860 4073
 Fax 209 368 4136
www.americanmastertech.com

StatLab
 2090 Commerce Drive
 McKinney, TX 75069 USA
 Tel 800 442 3573
 Fax 972 436 1369
www.statlab.com

1. A.F.I.P. Manual: pg 111. Strumia: Laboratory and Clinical Medicine; 930.
2. With modifications by AMTS R&D Department, 1979-2019.

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