Human Colon

Light Green Counterstain

100ml Item#: STLGC100

Liter Item#: STLGC1T

Pint Item#: STLGCPT

Gallon Item#: STLGC1GAL

Control Slide(s)
Item#

Candida Albicans CSC0525P

Pneumocystis CSP0125P

Required Components (Not Included)

Appropriate stain kit, such as GMS (Item#: KTGMS)

or PAS (Item#: KTPAS)

PRINCIPLE AND RESULTS: This stain is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to aid in the identification fungi and pneumocystis by staining background tissue green when it is used with a GMS or PAS stain procedure.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

G.M.S. STAIN KIT PROCEDURE (ITEM#: KTGMS):

METHENAMINE/SILVER SOLUTION PREPARATION: Prepare solution at time of use. Solution expires after one use. Methenamine/Silver Compound may soften during transit; this will not affect product performance. If desired, Methenamine/Silver Compound can be re-hardened by placing vials in a freezer. Methenamine/Silver Solution reacts with metals. Do not use metal forceps, racks, or tools with the solution.

#	Action	Amount	Chemical/Reagent	Details
1	Add	2ml	Distilled Water	Into a vial of Methenamine/Silver Compound. Break up chunks with an applicator stick. Pour contents of the vial into a chemically cleaned container or new/unused plasticware.
2	Add	46ml	Distilled Water	Mix thoroughly.
3	Add	2ml	Sodium Borate Concentrate	Mix for 2 minutes.

STAINING PROCEDURE: Preheating Required. See step 6 for more information. For microwave procedure, please see Item#: KTGMS instructions.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	Once complete, rinse in running DI water (1 minute) and continue.
4	Immerse	5% Chromic Acid	--	60	--	Once complete, rinse in running DI water (1 minute) and continue.
5	Immerse	1% Sodium Bisulfite	--	1	--	Immerse into preheated solution. Once complete, rinse in running DI water (1 minute) and continue.
6	Waterbath	Methenamine/Silver Solution	62°	15-25	--	Once complete, rinse in running DI water (1 minute) and continue.
7	Immerse	0.3% Gold Chloride	--	1	--	Once complete, rinse in running DI water (1 minute) and continue.
8	Immerse	2% Sodium Thiosulfate	--	1	--	Once complete, rinse in running DI water (1 minute) and continue.
9	Immerse	Light Green Counterstain	--	--	30-60	
10	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	Fresh changes. 1 minute each change.
11	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
12	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB; Theory and Practice of Histotechnology; 1980, 245 - 246.
2. Sale GE: Rapid Methenamine Silver Stain; Arch Pathol Lab Med: 102; 1978, 351 - 352.
3. With modifications by AMTS R&D Department, 1979-2019.

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		

CONTACT INFORMATION:

American MasterTech Scientific

1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab

2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. Sheehan DC Hrapchak BB; Theory and Practice of Histotechnology; 1980, 245 - 246.
2. Sale GE: Rapid Methenamine Silver Stain; Arch Pathol Lab Med: 102; 1978, 351 - 352.
3. With modifications by AMTS R&D Department, 1979-2019.

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