Human Spleen

Methylene Blue Stain Procedure

100ml Item#: STMBL100

Liter Item#: STMBLLT

Pint Item#: STMBLPT

Gallon Item#: STMBLGAL

Control Slide(s)	Item#
Acid Fast Bacteria	CSA0125P

Required Components (Not Included)
Advanced Carbol Fuchsin™
1% Acid Alcohol

PRINCIPLE AND RESULTS: This stain is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to aid in the identification of acid fast organisms. Methylene Blue Stain colors background tissue and nuclei blue and can be used as the counterstain for Kinyoun's acid fast stain and Ziehl-Neelsen stain procedures.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 6µm tissue section.

STORAGE AND USAGE NOTES: Filter Advanced Carbol Fuchsin™ when a shiny gold layer appears on the surface. Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

ACID FAST BACTERIA STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	Immerse	Advanced Carbol Fuchsin™	--	10-15	--	Once complete, rinse in running tap water (1 minute) and continue.
5	Immerse	1% Acid Alcohol	--	--	3-5	Differentiate with agitation until background is pale pink. Use a quick tap water rinse and microscope to check for correct differentiation. Once complete, rinse in running tap water (1 minute) and continue.
6	Immerse	Methylene Blue Stain	--	--	5-15	
7	Rinse	Running DI Water	--	--	3-5	
8	Dehydrate	Absolute Alcohol, 2 changes	--	--	15	15 seconds each change.
9	Clear	Xylene or Substitute	--	1	--	Or as required if using a xylene substitute.
10	Coverslip	Permanent Mounting Media	--	--	--	

1. A.F.I.P. Manual, pg. 205; Mallory Pathological Technique, pg. 275
2. With modifications by AMTS R&D Department, 1979-2019.

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		
 GHS07	Irritant	 GHS08	Health Hazard			

CONTACT INFORMATION:

American MasterTech Scientific
 1330 Thurman St.
 Lodi, CA 95240 USA
 Tel 800 860 4073
 Fax 209 368 4136
www.americanmastertech.com

StatLab
 2090 Commerce Drive
 McKinney, TX 75069 USA
 Tel 800 442 3573
 Fax 972 436 1369
www.statlab.com

1. A.F.I.P. Manual, pg. 205; Mallory Pathological Technique, pg. 275
2. With modifications by AMTS R&D Department, 1979-2019.

MULTILANGUAGE PROCEDURE

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COMPOSANTS NON INCLUS: Advanced Carbol Fuchsin™, 1% Acid Alcohol

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