Human Bone Marrow

## May Grunwald Stain (Stock) (Jenner Stain) Procedure

0:05

100ml Item#: N/A

Liter Item#: STMGRLT

Pint Item#: STMGRPT

Gallon Item#: STMGRGAL

### Control Slide(s)

### Item#

### Required Components (Not Included)

Bone Marrow CSB0325P

Phosphate Buffer pH 6.8

Helicobacter Pylori CSH0125P

Giemsa Stock

Mast Cell CSM0125P

 0.25% Acetic Acid (*Mast Cells only*)

Toxoplasma Gondii CST0325P

**PRINCIPLE AND RESULTS:** This stain is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify bone marrow cells, *Helicobacter pylori*, mast cells, rickettsias, and *Toxoplasma gondii*. Nuclei stain blue to violet, rickettsias intense reddish-purple, cytoplasm light blue, collagen and muscle pale pink, erythrocytes gray, yellow or pink, *Helicobacter pylori* blue, mast cells dark blue with red granules, and *Toxoplasma gondii* blue.

**SPECIMEN CRITERIA:** Appropriately fixed, paraffin-embedded, 5µm tissue section, or air dried smear.

**STORAGE AND USAGE NOTES:** Store/Use component according to the temperature and expiration on the label.

**PRECAUTIONS:** For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

**WORKING MAY-GRUNWALD STAIN PREPARATION:** Prepare solution at time of use. Solution expires after one use.

| # | Action | Amount | Chemical/Reagent        | Details                         |
|---|--------|--------|-------------------------|---------------------------------|
| 1 | Add    | 25ml   | May-Grunwald Stock      | Into a coplin jar or container. |
| 2 | Add    | 25ml   | Phosphate Buffer pH 6.8 | Mix thoroughly.                 |

**WORKING GIEMSA STAIN PREPARATION:** Prepare solution at time of use. Solution expires after one use.

| # | Action | Amount | Chemical/Reagent        | Details                         |
|---|--------|--------|-------------------------|---------------------------------|
| 1 | Add    | 50ml   | Phosphate Buffer pH 6.8 | Into a coplin jar or container. |
| 2 | Add    | 2.5ml  | Giemsa Stock            | Mix thoroughly.                 |

**STAINING PROCEDURE:** If staining a smear, begin procedure at step #3.

| # | Action                                                                                        | With                            | Heat °C | Time |      | Details                                                                 |
|---|-----------------------------------------------------------------------------------------------|---------------------------------|---------|------|------|-------------------------------------------------------------------------|
|   |                                                                                               |                                 |         | Mins | Secs |                                                                         |
| 1 | Deparaffinize                                                                                 | Xylene or Substitute, 2 changes | --      | 5    | --   | 5 minutes each change or as required if using a xylene substitute.      |
| 2 | Rinse                                                                                         | Absolute Alcohol, 3 changes     | --      | 1    | --   | 1 minute each change or as required if using graded alcohols.           |
| 3 | Rinse                                                                                         | Running DI Water                | --      | 1    | --   |                                                                         |
| 4 | Immerse                                                                                       | Working May-Grunwald Stain      | --      | 6    | --   | Or flood on slide rack. Agitate every few minutes for optimum staining. |
| 5 | Rinse                                                                                         | Phosphate Buffer pH 6.8         | --      | --   | 5-10 | Or flood on slide rack until no stain runs off slide.                   |
| 6 | Immerse                                                                                       | Working Giemsa Stain            | --      | 13   | --   | Or flood on slide rack. Agitate every few minutes for optimum staining. |
| 7 | Rinse                                                                                         | Phosphate Buffer pH 6.8         | --      | --   | 5-10 | Or flood on slide rack until no stain runs off slide.                   |
| 8 | To stain for mast cells, follow procedure (8a) below. For all other staining, skip to step 9. |                                 |         |      |      |                                                                         |

#### 8a. Mast Cell Procedure

|         |                   |    |    |   |                                                                                                                                                                                                  |
|---------|-------------------|----|----|---|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immerse | 0.25% Acetic Acid | -- | -- | 1 | 3-5 dips for 1 second each until desired background intensity is achieved. Use a quick Phosphate Buffer pH 6.8 rinse and microscope to check for correct differentiation and continue at step 9. |
|---------|-------------------|----|----|---|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|    |           |                               |    |      |    |                                  |
|----|-----------|-------------------------------|----|------|----|----------------------------------|
| 9  | Immerse   | Phosphate Buffer pH 6.8       | -- | 3    | -- | Or flood on slide rack.          |
| 10 | Air Dry   | On Benchtop, or in Empty Bath | -- | 5-10 | -- | Until section is completely dry. |
| 11 | Clear     | Xylene or Substitute          | -- | --   | 1  | 5-10 dips for 1 second each.     |
| 12 | Coverslip | Permanent Mounting Media      | -- | --   | -- |                                  |

1. A.F.I.P. Manual: pg 111. Strumia: Laboratory and Clinical Medicine; 930.
2. With modifications by AMTS R&D Department, 1979-2019.

**CE MARKINGS AND DESIGNATIONS:**

|                                                                                            |                                    |                                                                                            |                                   |                                                                                            |                |                                                                                                                  |
|--------------------------------------------------------------------------------------------|------------------------------------|--------------------------------------------------------------------------------------------|-----------------------------------|--------------------------------------------------------------------------------------------|----------------|------------------------------------------------------------------------------------------------------------------|
| <b>REF</b>                                                                                 | Catalogue Number                   |           | Temperature Limitation            |           | Manufacturer   | American MasterTech Scientific<br>1330 Thurman St.<br>Lodi, CA 95240 USA<br>Tel 800 860 4073<br>Fax 209 368 4136 |
| <b>LOT</b>                                                                                 | Batch Code                         |           | Use By                            | <b>EC REP</b>                                                                              | Representative | Emergo Europe<br>Prinsessegracht 20<br>2514 AP The Hague<br>The Netherlands                                      |
| <b>IVD</b>                                                                                 | In Vitro Diagnostic Medical Device |           | Consult Instructions Prior to Use | <b>CE</b>                                                                                  |                |                                                                                                                  |
| <br>GHS02 | Flammable                          | <br>GHS06 | Toxic                             | <br>GHS08 | Health Hazard  |                                                                                                                  |

**CONTACT INFORMATION:**

**American MasterTech Scientific**  
 1330 Thurman St.  
 Lodi, CA 95240 USA  
 Tel 800 860 4073  
 Fax 209 368 4136  
[www.americanmastertech.com](http://www.americanmastertech.com)

**StatLab**  
 2090 Commerce Drive  
 McKinney, TX 75069 USA  
 Tel 800 442 3573  
 Fax 972 436 1369  
[www.statlab.com](http://www.statlab.com)

1. A.F.I.P. Manual: pg 111. Strumia: Laboratory and Clinical Medicine; 930.
2. With modifications by AMTS R&D Department, 1979-2019.

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