Human Uterus

One-Step Trichrome Green & Red Stain Procedure

0:05

100ml Item#: STOSTG100

Liter Item#: STOSTGLT

Pint Item#: STOSTGPT

Gallon Item#: STOSTGGAL

Control Slide(s)	Item#
Heart	CSH0725P
Skeletal Muscle	CSS0725P
Trichrome	CST0125P
Uterus	CSU0325P

Required Components (Not Included)
 Bouin's Fluid or Bouin's 2000™
 Modified Mayer's Hematoxylin


PRINCIPLE AND RESULTS: This stain is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify collagen and muscle. Cytoplasm, keratin, muscle, intercellular fibers are stained red, collagen and mucus green, and nuclei blue.

SPECIMEN CRITERIA AND CONTROL TISSUE: Appropriately fixed, paraffin-embedded, 4-6µm tissue section.

STORAGE AND USAGE NOTES: One Step Trichrome Stain may require filtering. Store/Use component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

STAINING PROCEDURE: Preheating Required. See step 4 for more information.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	Waterbath	Bouin's Fluid or Bouin's 2000™	56°	60	--	Mordant in preheated Bouin's Fluid for 60 minutes. For Bouin's 2000™, mordant overnight at room temperature. Once complete, rinse in running tap water (3 minutes) or until tissue is colorless and continue.
5	Immerse	Modified Mayer's Hematoxylin	--	7	--	Once complete, rinse in running tap water (3 minutes) and continue.
6	Immerse	One Step Trichrome Stain (Green & Red)	--	3-8	--	Once complete, rinse in running tap water (5 seconds) and continue.
7	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
8	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
9	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 191.
2. With modifications by AMTS R&D Department, 1979-2019.

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		

CONTACT INFORMATION:
American MasterTech Scientific

1330 Thurman St.
 Lodi, CA 95240 USA
 Tel 800 860 4073
 Fax 209 368 4136
www.americanmastertech.com

StatLab

2090 Commerce Drive
 McKinney, TX 75069 USA
 Tel 800 442 3573
 Fax 972 436 1369
www.statlab.com

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 191.
2. With modifications by AMTS R&D Department, 1979-2019.

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