Human Breast Tissue

Picric Acid / Acetone (Traditional Formulation) Stain Procedure

100ml Item #: N/A

Liter Item#: N/A

Pint Item #: STPAAPT

Gallon Item#: STPAAGAL

Control Slide(s)	Item#
Gram +/- Gram -	CSG0225P

Required Components (Not Included)	
Gentian Violet	1% Safranin O
Universal Iodine Solution™	Gram's Decolorizer

PRINCIPLE AND RESULTS: This stain is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens or smears (in vitro) to assist in the identification of gram positive and gram negative bacteria. Picric Acid Acetone differentiates tissue colors by removing excess Safranin O stain while staining the background yellow.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	Immerse	Gentian Violet	--	2	--	Or flood slide on slide rack. Once complete, rinse in running tap water (1 minute) and continue.
5	Immerse	Universal Iodine Solution™	--	1	--	Or flood on slide rack. Once complete, rinse in running tap water (1 minute) and continue.
6	Decolorize	Gram's Decolorizer	--	--	5-10	Until color does not bleed from tissue. Once complete, rinse in running tap water (5-10 seconds) and continue.
7	Immerse	1% Safranin O	--	3	--	Once complete, rinse in running tap water (5-10 seconds) and continue.
8	Immerse	Picric Acid/Acetone	--	--	1	1-5 dips for 1 second each, until correct differentiation is achieved.
9	Dehydrate	Absolute Alcohol, 2 changes	--	--	1-2	Fresh changes. Dip for 1-2 second each change.
10	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
11	Coverslip	Permanent Mounting Media	--	--	--	

1. Carson FL Hladik C: Histotechnology, A Self-Instructional Text; 3rd Edition, 2009, 231-232.
2. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 234-235.
3. With modifications by AMTS R&D Department, 1979-2019.

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		
 GHS02	Flammable	 GHS07	Irritant			

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. Carson FL Hladik C: Histotechnology, A Self-Instructional Text; 3rd Edition, 2009, 231-232.
2. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 234-235.
3. With modifications by AMTS R&D Department, 1979-2019.

MULTILANGUAGE PROCEDURE

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COMPOSANTS INCLUS: Gentian Violet, Universal Iodine Solution™, Gram's Decolorizer, 1% Safranin O

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