Human Uterus

P.T.A.H. (Phosphotungstic Acid Hematoxylin) Stain Procedure

100ml Item#: STPTA100

Liter Item#: N/A

Pint Item#: STPTAPT

Gallon Item#: N/A

Control Slide(s)
Item#
Required Components (Not Included)

Heart CSH0725P

10% Zinc Chloride

Skeletal Muscle CSS0725P

5% Ferric Ammonium Sulfate

Trichrome CST0125P

Uterus CSU0325P


PRINCIPLE AND RESULTS: This stain is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify collagen, fibrin and striated muscle without using Zenker's Fixative with Mercuric Chloride as a mordant. Fibrin and muscle stain blue, collagen brownish-red to red, and nuclei blue.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

STORAGE AND USAGE NOTES: Store/Use component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

STAINING PROCEDURE: Preheating Required. See step 6 for more information. For a microwave procedure, please see instructions for item#: KTPTA.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohols, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse 	Running DI Water	--	1	--	
4	Immerse	10% Zinc Chloride Solution	--	180	--	Once complete, rinse in running DI water (1 minute) and continue.
5	Immerse	5% Ferric Ammonium Sulfate Solution	--	10	--	Without rinsing, continue to next step.
6	Waterbath	P.T.A.H. Stain	65°	120	--	Immerse into preheated solution. For room temperature staining without a waterbath, increase time to 12-24 hours while periodically agitating.
7	Differentiate	95% Reagents Alcohol	--	--	30-60	Check using microscope for correct differentiation.
8	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
9	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
10	Coverslip	Permanent Mounting Media	--	--	--	

1. A.F.I.P. Laboratory Methods in Histotechnology: 1992; 98 - 99.
2. With modifications by AMTS R&D Department, 1979-2019.

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		

CONTACT INFORMATION:

American MasterTech Scientific

1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab

2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. A.F.I.P. Laboratory Methods in Histotechnology: 1992; 98 - 99.
2. With modifications by AMTS R&D Department, 1979-2019.

MULTILANGUAGE PROCEDURE

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COMPOSANTS NON INCLUS: 10% Zinc Chloride Solution, 5% Ferric Ammonium Sulfate Solution

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