

Instructions for Use

Acid Fast Bacteria (AFB) Control Slides

CATALOG NUMBER	DESCRIPTION	UNIT OF MEASUREMENT
CS003-25	AFB Control Slides, Non-adhesion, 1 stained/24 unstained	1 set
CS003-PS-25	AFB Control Slides, Adhesion, 1 stained/24 unstained	1 set
CS003-100	AFB Control Slides, Non-adhesion, 1 stained/99 unstained	1 set
CS003-PS-100	AFB Control Slides, Adhesion, 1 stained/99 unstained	1 set
CS101-PS-25	AFB +/- Control Slides, Adhesion, 1 stained/24 unstained	1 set
CS101-PS-100	AFB +/- Control Slides, Adhesion, 1 stained/99 unstained	1 set

INTENDED USE

Qualitative control for acid fast bacteria.

Please contact productsupport@statlab.com with any additional questions.

TISSUE


Known Positive Tissue

FIXATION

10% Buffered Formalin

TECHNICAL INFORMATION

Paraffin sections cut at 5 microns



USE WITH MASTERTECH A.F.B. STAIN KIT

Kit Size	Part Number
100 mL	KTAFB
Pint	KTAFBPT
Liter	KTAFBLT

MASTERTECH A.F.B. STAIN KIT PROCEDURE

1. Deparaffinize and hydrate to running tap water.
2. Immerse in Advanced Carbol Fuchsin for 10-15 minutes.
3. Rinse in running tap water for 1 minute.
4. Dip in 1% Acid Alcohol for 3-5 seconds with agitation until background is pale pink. Use a quick tap water rinse and microscope to check for correct differentiation.
5. Rinse in running tap water for 1 minute.
6. Counterstain in Light Green Counterstain for 1 minute. DO NOT RINSE IN WATER.
7. Dehydrate in 3 changes of Absolute Alcohol, 1 minute each change.
8. Clear in 3 changes of Xylene, 1 minute each change.
9. Coverslip with Acrymount Plus or any other acceptable mounting medium.

RESULTS

All acid fast bacteria – bright red

Other tissue – light green

This protocol is a recommendation only. Modifications may be necessary.



StatLab Medical Products
2090 Commerce Drive, McKinney, TX 75069
StatLab.com Phone: +1 800-442-3573

