

**GOMORI'S METHOD FOR IRON**

<b>PURPOSE:</b>	For In Vitro Diagnostic Use: Intended for the qualitative demonstration of ferric iron in tissue.
<b>PRINCIPLE:</b>	Loosely bound ferric iron in tissue will readily bond with Potassium Ferrocyanide in an acid medium forming an insoluble blue pigment (Prussian blue).
<b>CONTROL:</b>	Normal liver <i>Control Slides can be purchased from Histology Control Systems. See inside back cover, Item# cs007.</i>
<b>SPECIMEN PREPARATION:</b>	Formalin fixed, paraffin embedded sections cut at 6 micrometers
<b>SOLUTIONS:</b>	1. Hydrochloric Acid 20% Aqueous Item# s1832 2. Potassium Ferrocyanide 10% Item# s1873 3. Nuclear Fast Red Kernechtrot 0.1% Item# s248  <i>Solutions can be purchased separately from Poly Scientific.</i>
<b>NOTES:</b>	Must be prepared fresh each time. <u>Hydrochloric Acid-Potassium Ferrocyanide Working Solution:</u> Hydrochloric Acid 20%.....50 mL Potassium Ferrocyanide 10%.....50 mL
<b>REFERENCE:</b>	Gomori, G. "Gomori's Method for Iron". <u>Amer J. Path.</u> 12:655-663. 1936.

**STAINING PROCEDURE:**

1. Deparaffinize and hydrate to distilled water.
2. Place in Hydrochloric Acid-Potassium Ferrocyanide Working Solution for 15 minutes.
3. Rinse in distilled water.
4. Counterstain in Nuclear Fast Red Kernechtrot 0.1% for 2 minutes.
5. Rinse in distilled water.
6. Alcohol 95%, 2 changes.
7. Alcohol 100%, 2 changes.
8. Xylene, 2 changes.
9. Mount with Poly Mount (Item# s2153) or any other acceptable mounting medium.

**RESULTS:**

Iron Pigments.....	Bright Blue
Cytoplasm .....	Light Pink
Nuclei.....	Red

*Poly Scientific R&D Corp.*

Revision: B-18

