

**GIEMSA RAPID STAIN METHOD FOR HELICOBACTER PYLORI**

<b>PURPOSE:</b>	For In Vitro Diagnostic Use: Intended for the qualitative demonstration of Helicobacter Pylori.
<b>PRINCIPLE:</b>	The Giemsa Stain is a combination of acid and basic dyes. Differential staining is attributed to relative charge of cells, dye size and the pH of the solutions.
<b>CONTROL:</b>	Tissue known positive for Helicobacter Pylori <i>Control Slides can be purchased from Histology Control Systems. See inside back cover, Item# cs033.</i>
<b>SPECIMEN PREPARATION:</b>	Formalin fixed, paraffin embedded sections cut at 4-5 micrometers
<b>SOLUTIONS:</b>	<p>1. Giemsa Stain Stock Solution Item# s195 <u>Giemsa Working Solution:</u> Prepare fresh, filter before use, do not reuse.</p> <p style="margin-left: 20px;">a. Giemsa Stock Solution.....10 mL b. Distilled Water.....20 mL</p> <p>2. Wright Stain Buffer pH 6.4 Giordano Item# s2123 <i>Solutions can be purchased separately from Poly Scientific.</i></p>
<b>NOTES:</b>	
<b>REFERENCE:</b>	Luna, Lee G. <u>Histopathologic Methods and Color Atlas of Special Stains and Tissue Artifacts</u> . American Histolabs Inc. Gaithersburg, MD. 1992. p. 216.

**STAINING PROCEDURE:**

1. Deparaffinize and hydrate to distilled water.
2. Place slides in Giemsa Stain Working Solution for 3 minutes.
3. Quickly dip slides twice in Wright Stain Buffer pH 6.4 Giordano.
4. 3 changes of Absolute Alcohol.
5. 3 changes of Xylene.
6. Mount with Poly Mount (Item# s2153) or any other acceptable mounting medium.

**RESULTS:**

Helicobacter Pylori..... Blue  
Background ..... Various Shades of Blue

*Poly Scientific R&D Corp.*

Revision: B-18

