Animal Liver

Chandler's Precision® Reticulum Stain Kit Procedure

100ml Kit Item #: KTCPR

Liter Kit Item#: KTCPRLT

Pint Kit Item #: KTCPRPT

Gallon Kit Item#: N/A

Control Slide(s)	Item#
Reticulum	CSR0125P

Included Components

Type 1 Silver Solution	20% Formalin
3% Sodium Hydroxide	0.1% Gold Chloride
1% Potassium Permanganate	5% Sodium Thiosulfate
3% Potassium Metabisulfite	Nuclear Fast Red Stain
3% Ferric Ammonium Sulfate	

Components Not Included:

Concentrated Ammonium Hydroxide (Item#:05870-16)

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify reticulum fibers. Reticulum fibers stain gray to black and nuclei pink to red.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 5µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

AMMONIACAL SILVER SOLUTION PREPARATION: Prepare solution at time of use. Solution expires after one use. Ammoniacal Silver Solution reacts with metals. Do not use metal forceps, racks, or tools with the solution.

#	Action	Amount	Chemical/Reagent	Details
1	Add	1 vial	Type 1 Silver Solution	Into a chemically cleaned container or new/unused plasticware.
2	Add	17ml	Distilled Water	Mix thoroughly.
3	Add	1ml	Concentrated Ammonium Hydroxide (Not Included)	Drop by Drop. Mix thoroughly; solution will turn brown, then colorless.
4	Add	10ml	3% Sodium Hydroxide	Mix thoroughly. If no precipitate remains, skip to Step 6.
5	Add	Drops	Concentrated Ammonium Hydroxide (Not Included)	Drop by Drop while mixing, until no precipitate remains.
6	Add	30ml	Distilled Water	Mix thoroughly.

STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Immerse	1% Potassium Permanganate	--	5-10	--	Once complete, rinse in running DI water (1 minute) and continue.
5	Immerse	3% Potassium Metabisulfite	--	--	30-60	Until section is colorless. Once complete, rinse in running DI water (1 minute) and continue.
6	Immerse	3% Ferric Ammonium Sulfate	--	10	--	Once complete, rinse in running DI water (1 minute) and continue.
7	Immerse	Ammoniacal Silver Solution	--	2-3	--	Once complete, rinse in running DI water (1 minute) and continue.
8	Immerse	20% Formalin	--	1	--	Once complete, rinse in running DI water (1 minute) and continue.
9	Immerse	0.1% Gold Chloride	--	3-5	--	Once complete, rinse in running DI water (1 minute) and continue.
10	Immerse	5% Sodium Thiosulfate	--	1-2	--	Once complete, rinse in running DI water (1 minute) and continue.
11	Immerse	Nuclear Fast Red Stain	--	3-5	--	Once complete, rinse in running DI water (1 minute) and continue.
12	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
13	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
14	Coverslip	Permanent Mounting Media	--	--	--	

- Gomori: A Modification of the Silver Impregnation Method of Staining Reticular Fibers; American Journal of Clinical Pathology, 21:1951, 897-899
- With modifications by AMTS R&D Department, 1979-2018.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. Standard Ammoniacal Silver Solution preparation procedure yields approximately 61ml of solution and must be scaled up to accommodate desired autostainer bath size. A minimum of 9 baths are required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 20 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	60	240	4	250ml Glass Stain Dish	4	120	30
30ml Glass Coplin Jar	40	320	8	200ml Bath Autostainer	6	180	30
40ml Hellendahl Jar	30	480	16	400ml Bath Autostainer	3	150	50

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	StatLab Medical Products 2090 Commerce Drive McKinney, TX 75069 www.statlab.com 800-442-3573
LOT	Batch Code		Use By			
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	 GHS05	Corrosive	
 GHS07	Irritant	 GHS08	Health Hazard	 GHS09	Environmentally Damaging	

- Gomori: A Modification of the Silver Impregnation Method of Staining Reticular Fibers; American Journal of Clinical Pathology, 21:1951, 897-899
- With modifications by AMTS R&D Department, 1979-2018.

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