DIAPATH

Constant Innovation in Anatomic Pathology

INSTRUCTION FOR USE

CPD06

EN

IT

FR

ES

DE

CytoPath® Disc 6 mm

Device for cell block preparation

EN

Manufacturer

Diapath S.p.A.

Via Savoldini, 71 - 24057 Martinengo (BG) - Italy

Tel. 0363 986411

Fax 0363 948000

www.diapath.com

customerservice@diapath.com

EMDN code: W0202059085

UDI code: 08033895010191

Code	Test	Components
CPD06	20	6 mm Disc matrices

Intended use

Single-use, non-sterile in vitro diagnostic device for clinical use, for the collection, uptake and preparation of cell block material *in vitro*, previously obtained in accordance with the method of reference for the type of histological sample, which can then be used in the various processing stages envisaged in an anatomical pathology laboratory, including molecular biology and immunohistochemistry techniques. The device is compatible with commercially available cytology fixatives.

CytoPath®- Disc 6mm is packaged, ready for use.

The IVD CPD06 device must only be used by experienced and qualified personnel, by means of manual operation, free of toxic reagents. The device is a circular disc of polymeric material used in conjunction with consumables, such as falcon, eppendorf, pipettes, and instruments, such as centrifuge, processor, embedding station, microtome, stainer.

The CytoPath®- Disc 6mm IVD device allows the uptake of the cytological pellet, once it has been inserted into the falcon, from which the supernatant has already been separated and removed, in accordance with the laboratory's procedures, making it available for subsequent anatomical pathology processing and enabling the preparation of a cell block.

The device does not come into contact with the patient as it is a diagnostic aid; it is used in the sample preparation stage for research or diagnostic purposes.

The device does not provide specific results as it is involved in the treatment process in order to prepare the sample for the next step; therefore, the device is an aid to diagnosis and does not provide any interpretable analytical results.

For further information, warnings and precautions for use, please see the website.

Description of the reagents and any limitations of use

Device	Stability of closed product	Ready for use	Storage temperature
CPD06	2 years	YES	15°C -30°C

Information relating to interfering substances that may affect the performance of the device

Not applicable

List of materials provided and specific materials required but not provided

The device is packaged, ready to use, in units of 6 mm circular Disc matrices.

The following materials are required for the analysis, but are **not supplied**:

- Consumables and instruments required for cytological sampling and for the centrifugation stage (centrifuge, pipettes, falcon, eppendorf, etc.)
- Loops, pipettes or tweezers to retrieve the 6 mm Disc matrix once the cytological material has been absorbed
- Consumables and tools required for cell block processing and the subsequent steps of embedding, cutting, staining and mounting
- All necessary instrumentation and kits for subsequent molecular biology, immunohistochemistry and immunocytochemistry analyses

Information on the combined use of the device

The device is used in conjunction with:

- Instruments and consumables required for cell block processing and the subsequent steps of embedding, cutting, staining and mounting
- All necessary instrumentation for carrying out additional molecular biology, immunohistochemistry and immunocytochemistry analyses

The device may only be used in combination with the devices listed above. Any use other than intended is only permitted after internal laboratory validation.

Functional characteristics

- Type of preparation: Cytological samples
- Ready-to-use device; no preparation prior to uptake of the cell pellet is required
- Indicative time for preparation of cytological material to be processed: approx. 15 minutes for absorption into the Disc 6 mm matrix
Requirements: The absorption time of the cytological material within the Disc matrix must be no less than 10 minutes, followed by fixation in buffered formalin 10% v/v - 4% w/v or other fixative for at least 4 hours before going ahead with the processing.

Test principle

CytoPath® Disc 6 mm makes it possible to obtain cell blocks from gynaecological and non-gynaecological cytological material previously obtained in accordance with the prescribed cytological reference methods.

The operation of CytoPath® Disc 6 mm is based on the absorption of the cell pellet. The polymer matrix that has absorbed the cytological pellet will then be processed in accordance with histological preparation methods.

Instructions for use

METHOD 1

- 1) The cytological sample taken is placed inside a falcon or eppendorf and centrifuged at 600 rev for 10 minutes by means of a centrifuge, the supernatant is removed with a pipette or other laboratory method, leaving a deposit of cells (pellet) at the base of the container (Fig 1).

Figure 1. Within a falcon it is necessary to separate the cytological pellet from the supernatant by centrifugation

- 2) Insert the CytoPath® Disc 6 mm device inside the falcon or eppendorf, only after removing the supernatant, so that it settles on the conical bottom of the container in order to come into direct contact with the cell pellet (Fig.2).

Figure 2. Remove the supernatant, insert the device, wait 10 minutes for absorption of cellular material into the disc

- 3) Close the falcon or eppendorf inside which the device is located and let it absorb for 10 min.
- 4) Remove the CytoPath® Disc 6 mm from the falcon with a plastic pasteur pipette or similar device, and place it inside filter bag or filter paper (this minimises the loss of cellular material in subsequent steps). Insert the filter bag or filter paper inside a histology cassette
- 5) Fix in buffered formalin 10% v/v - 4% w/v for at least 4 hours, as inadequate fixation can sometimes make subsequent cutting of sections difficult and affect the quality of the immunohistochemistry (Fig.3)


Figure 3. Once the cassette is closed, place it inside a container with buffered formalin 10% v/v - 4% w/v for at least 4 hours.

- 6) Proceed with the processing step using a tissue processor (Fig.4)

Figure 4. Go ahead with processing with the tissue processor

- 7) Pull out CytoPath® Disc 6 mm and include the device by orienting it in such a way that the part with the most cells faces the lower surface of the mould; proceed with the microtome cutting, staining and mounting steps or other method prescribed by the laboratory (Fig. 5)

NOTE: For fixation and staining requirements please refer to the target laboratory routine; therefore, users must evaluate their own fixation and staining protocols for an optimal result.

Figure 5. Include the block, cut, stain and assemble the sections.

METHOD 2

In the case of low-cell and/or low-volume samples, the following method is recommended in order to achieve sample enrichment.

- 1) Collect the sample by means of a pasteur pipette, micropipette or syringe, so that it can be aspirated and can be easily expelled.
- 2) Place the CytoPath® Disc 6 mm on a clean surface, (e.g. use a clean support slide)
- 3) Deposit the sample directly onto the synthetic matrix, depositing it slowly so that the matrix can absorb all the material. (Fig.6)


Figure 6. Sample drop deposited directly onto CytoPath® Disc 06 mm

- 4) Insert CytoPath® Disc 6 mm inside a filter bag or filter paper (this minimises the loss of cellular material in subsequent steps). Insert the filter bag or filter paper inside a histology cassette
- 5) Fix in buffered formalin 10% v/v - 4% w/v for at least 4 hours, as inadequate fixation can sometimes make subsequent cutting of sections difficult and affect the quality of the immunohistochemistry (Fig.3)
- 6) Proceed with the processing step using a tissue processor (Fig.4)
- 7) Pull out CytoPath® Disc 6 mm and include the device by orienting it in such a way that the part with the most cells faces the lower surface of the mould; proceed with the microtome cutting, staining and mounting steps or other method prescribed by the laboratory (Fig. 5)

Please refer to the Bibliography section

It is neither a self-diagnostic device nor is it intended for decentralised analysis. This device may only be used by qualified and suitably trained personnel, since it is mainly intended for health professionals and operators working in specialised medical facilities. Occupational safety guidelines must be applied in accordance with the applicable regulations. The instruments used for making the diagnosis must be suitable for diagnostic use in the laboratory. A diagnosis must only be made by authorised, qualified and trained personnel.

Analytical performance characteristics

Not applicable

Conditions for sampling, handling and preparation of the sample

The subsequent activities for preparation of the sample follow the normal workflow of a pathology laboratory.

Do not handle until you have read and understood all warnings.

Storage

Optimum storage conditions are at a temperature of 15 - 30°C and relative humidity of 20 - 80%.

Do not expose the device to environmental conditions other than those prescribed, as failure to comply with storage conditions would lead to a change in the chemical and physical characteristics of the device and, consequently, a non-compliant result.

If properly stored and in unopened packaging, the product is stable as stated in the section *Description of the reagents and any limitations of use*. Do not use after the expiry date.

If the product is not stored as instructed, its yield may vary and it will have to be tested by the operator.

Stability when handling the device

The use of the device is under the absolute control of the skilled and qualified operator.

The environmental conditions in the laboratory must be controlled:

- Temperature: +15°C / +30°C
- Relative humidity ≤80%

Disposal instructions / Safe final disposal of the device

For specific instructions on how to safely dispose of the device after use, please refer to the local regulations of the country in which the device is used.

It is good practice to dispose of the device as special waste in dedicated containers and in accordance with internal disposal procedures.

For information on the content of the reagents, please refer to the safety data sheet. If the product has expired and/or is unused, it must be disposed of in accordance with local waste regulations, based on the hazard classification on the label, after checking for any contamination. In some cases, analytical evaluation may be required so as to correctly classify the waste and determine the hazard characteristics to be attributed to it

Warnings, precautions and restrictions on the use of the device

The device is to be used only for its intended use.

For proper handling, it is recommended to:

- Wear personal protective equipment, such as a lab coat, gloves and goggles. Protect your eyes and face
- In order to obtain the expected result, it is advisable to follow the indications given in the protocol of use and the validated components within the kit
- Maintain a clean work environment and follow the cleaning protocols in the laboratory of use

Device malfunction or degradation

The device retains its properties throughout its estimated useful life if the storage requirements are guaranteed and complied with.

The device does not generate risks from any debris due to wear and tear, degradation products and processing residues.

In case of malfunction, isolate the device, replace it with a new one and dispose of the non-functioning kit. Should this occur, please contact the manufacturer by sending an e-mail to customerservice@diapath.com

Error	Possible cause	Possible corrective action
Fraudulent use beyond the expiry date	Inadequate device performance	Including appropriate warnings and markings on the label
Inadequate storage of the device		Inserting the appropriate symbols and storage warnings on the label

Warnings from the user, with the need to report any serious accident that occurred in relation to the device

In case of a serious incident, the event involving the device in question should be reported to the manufacturer or to the competent authorities, so that the evaluation of such accidents can be more easily coordinated and appropriate safety corrective actions can be taken.

Key of the symbols used in the labelling and instructions for use

Symbol	Symbol name	Explanatory notes
	Manufacturer	Indicates the manufacturer of the device
	Expiration date	Indicates the date after which the device must not be used
	Lot number	Indicates the manufacturer's batch necessary for batch identification
	Catalogue code	Indicates the code assigned by the manufacturer allowing its identification
	Distributor	Indicates the company distributing the medical device
	Not sterile	Indicates a device that has not undergone a sterilization process
	Keep away from light	Indicates a device that needs protection from light sources.
	Temperature limit	Indicates the temperature limits to which the device can be safely exposed

	Not reusable	Indicates a device that is to be used once, for single use only
	Refer to instructions	Indicates the need for the user to consult the instructions
	Attention	Indicates the need for the user to be cautious when operating the device or that the situation requires greater awareness to avoid undesirable consequences
	In vitro diagnostic medical device	Indicates a device intended to be used as an in vitro diagnostic medical device
	CE Marking	The CE mark certifies the conformity of the product to the applicable directives of the European Union
	Content sufficient for <n> tests	Indicates the number of tests that can be performed with the device
	Manufacturing date	Indicates the date of manufacture of the medical device
	Unique identifier of the device	This indicates a carrier containing the fixed and variable information related to the UDI

Bibliography

- Taccogna, Silvia, et al. "The Usefulness of CytoFoam-Core in the Diagnosis of Thyroid Metastasis From Colon Cancer: A Case Report." *Multidisciplinary Cancer Investigation* 4.4 (2020): 1-5.
- Cassano, Enrico, and Chiara Trentin. "Integrated breast biopsy for best radiological diagnosis of breast cancer." *Breast Cancer: Innovations in Research and Management* (2017): 317-331.
- Mayall, Frederick George, et al. "Cytology and cell-block immunohistochemistry of circulating tumour cells." *Cytopathology* 30.6 (2019): 620-627.
- Stevenson, Tracey, Manish Powari, and Christopher Bowles. "Evolution of a rapid onsite evaluation (ROSE) service for endobronchial ultrasound guided (EBUS) fine needle aspiration (FNA) cytology in a UK Hospital: A 7-year audit." *Diagnostic Cytopathology* 46.8 (2018): 656-662.

Revision table

Rev.	Date	Purpose of the revision
1st issue	18/03/2024	First issue
01	24/07/2025	Changes to the usage protocol
02	03/06/2026	Update to the 'Functional Features' section – amendment to the user guide

(*) If the request is due to a change, please provide a detailed description of the reason for revising the document and the amendments that have been made since the last revision; if available, please enclose recalls, non-conformities, NC AC AP statistics, etc.

You can find this document in the reserved area of the www.diapath.com website. If you want to download the document and keep up to date on the latest available version, you will need to log in after registering.

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Bibliographie

- Taccogna, Silvia, et al. « The Usefulness of CytoFoam-Core in the Diagnosis of Thyroid Metastasis From Colon Cancer : A Case Report. » *Multidisciplinary Cancer Investigation 4.4 (2020)* : 1-5.
- Cassano, Enrico et Chiara Trentin. « Integrated breast biopsy for best radiological diagnosis of breast cancer. » *Breast Cancer: Innovations in Research and Management (2017)* : 317-331.
- Mayall, Frederick George, et al. « Cytology and cell-block immunohistochemistry of circulating tumour cells. » *Cytopathology 30.6 (2019)* : 620-627.
- Stevenson, Tracey, Manish Powari et Christopher Bowles. « Evolution of a rapid onsite evaluation (ROSE) service for endobronchial ultrasound guided (EBUS) fine needle aspiration (FNA) cytology in a UK Hospital : A 7-year audit. » *Diagnostic Cytopathology 46.8 (2018)* : 656-662.

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