

This procedure is available in kit form from:  
**Poly Scientific R&D Corp.**  
70 Cleveland Avenue • Bay Shore, NY 11706  
(631) 586-0400

cat# cy020

**CRESYL VIOLET PROCEDURE (MOORE & BARR)**

**SOLUTIONS:** 1. Cresyl Violet 1% Aqueous

**STAINING PROCEDURE:**

1. Place in water for 5 minutes.
2. Stain in Cresyl Violet 1% Aqueous for 7 minutes.
3. Dehydrate in 95% Alcohol, 2 changes for 5 minutes (total).
4. Place in Absolute Alcohol for 5 minutes.
5. Clear in Xylene, 15 dips.
6. Mount with Poly Mount (cat# s 2153) or any other permanent mounting medium.

**RESULTS:**

Nuclei ..... Pale Pink  
Sex Chromatin ..... Deep Pink to Violet