

Sterile Slides (frosted, white yellow or plain ends)

Manufacturer/Supplier:

Scientific Device Laboratory
411 Jarvis Ave, Des Plaines, IL 60018 USA
General and Technical Information
Phone Number: 847-803-9495
Website: www.scientificdevice.com

Intended Use:

Sterile Slides are frosted, white, yellow or plain –ended slides measuring 1 inch x 3 inches. The slides are individually packaged and sealed in pouches and have been autoclaved. The sterile slides can be used in aseptic environments.

Summary and Explanation:

The slides have been tested for sterility with standard sterility techniques.

Components:

Sterile Slides are available in boxes of 50 individually bagged slides.

Warnings:

Do not use as sterile slide if pouch is ripped or compromised in any way.

Storage:

Store at room temperature (15°C to 30°)

Quality Control:

These Sterile Slides are quality controlled at Scientific Device Laboratory prior to shipping.

Limitations:

Do not use as sterile slide if pouch is ripped or compromised in any way.









Safety:

For *In vitro* diagnostic use.

Related Product Information:

Catalog # 050-S: Sterile Slides – Frosted End (50 slides per box)
Catalog # 050-1S: Sterile Slides – White End (50 slides per box)
Catalog # 050 – 2S: Sterile Slides – Yellow End (50 slides per box)
Catalog #050-SN: Sterile Slides – Plain End (50 slides per box)

SYMBOL LEGEND

	Catalog Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage Temp.)
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Quantity Included
	Dimensions

Supported Products Page **(QUALITY FILE ONLY)**

Revision History

CR NUMBER	REVISION
0410-002	00
1017-003	01
0819-001	02