

- Poly-l-lysine Coated Slides
- Silane Coated Slides
- Hybridization Coated Slides









Adhesive Microscope Slides

- Poly-l-lysine Coated Slides
- Silane Coated Slides
- Hybridization Coated Slides

Manufacturer/Supplier:

Scientific Device Laboratory
411 Jarvis Ave, Des Plaines, IL 60018 USA
General and Technical Information
Phone Number: 847-803-9495
Website: www.scientificdevice.com

SYMBOL LEGEND

	Catalog Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage Temp.)
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Quantity Included
	Dimensions

Intended Use:

This product is to be used with either tissue or smears that need an adhesive to adhere to the slide.

Summary and Explanation:

These adhesives have charged amino groups that adhere to both the glass surface and the sample being tested.

Components: (one of the following)

Poly-l-lysine
Silane (3-Aminopropyl) triethoxy-silane
Hybridization (*proprietary ingredients*)

Storage:

Store in a cool dry room.

Procedure:

Refer to the protocol that the laboratory has established.

Expected Results:

More than 95% of tissue or sample should adhere to microscope glass slide.

Quality Control:

These adhesive slides are quality controlled at Scientific Device Laboratory prior to shipping.

Limitations:

Keep away from wet warm areas.

Safety:

For *In vitro* diagnostic use. See MSDS for additional information.

Related Product Information:

Catalog # 067: Poly-l-lysine Coated Slides (minimum order of 5 gross)
Catalog # 068: Silane Coated Slides (minimum order of 5 gross)
Catalog # 063: Hybridization Coated Slides (minimum order of 5 gross).

- Poly-l-lysine Coated Slides
- Silane Coated Slides
- Hybridization Coated Slides

Supported Products Page **(QUALITY FILE ONLY)**

This Technical Insert is used for the following products:

SDL Prod ID	Description		
57	067 .	Poly-L-Lysine Adhesive Coated Slides	per gross (144)
	067 .-5	Poly-L-Lysine Adhesive Coated Slides-FISHER packaging	per case 5 gross (720)
59	068 .-5	Silane Adhesive Coated Slides-FISHER packaging	per case 5 gross (720)
	068 .	Silane Adhesive Coated Slides	per gross (144)
55	063	Hybredization Slide	per gross (144)

Revision History

CR NUMBER	REVISION
0410-002	00
1017-003	01