

Shrink Seals

Manufacturer/Supplier:

Scientific Device Laboratory
411 Jarvis Ave, Des Plaines, IL 60018 USA
General and Technical Information
Phone Number: 847-803-9495
Website: www.scientificdevice.com

Intended Use:

Shrink Seals are colorless, transparent cellulose bands used for sealing Petri dishes or test tubes.

Summary and Explanation:

Shrink Seals are used on plates to avoid drying of media. They are also good for safe and secure sealing of plates containing fungal, TB, or bacterial cultures that require long incubations.

Components:

Shrink Seals come in containers containing a preservative liquid. Shrink Seals are available in a variety of sizes and amounts.

Warnings:

Keep Shrink Seals wet in original container to prevent drying out. Make sure the top of the Shrink Seal container is securely closed.
Do not freeze.

Storage:

Store in original container with liquid in a cool place to prevent drying out of the seals.

Instructions for Use:

Remove Shrink Seals from container, lightly blot with paper towel. While Shrink Seals are still slightly wet, place seal around edge of Petri dish or test tube. Let edges of Shrink Seal overhang equally on both sides. Allow to air dry. Reseal container making sure lid is replaced tightly.

Safety:

For *In vitro* diagnostic use. See SDS for additional information.

Related Product Information:

Catalog # 261: Shrink Seals (92 x 25 mm) 400 seals per package
Catalog # 262: Shrink Seals (92 x 25 mm) 6,000 seals per package
Catalog # 269: Shrink Seals for Test Tubes (21 mm x 20 mm) 1,000 seals per package

SYMBOL LEGEND



Catalog Number



In Vitro Diagnostic Medical Device



Consult Instructions for Use (IFU)



Temperature Limitations (Storage Temp.)



Batch Code (Lot Number)



Use By (Expiration Date)



Quantity Included



Dimensions

Supported Products Page **(QUALITY FILE ONLY)**

This Technical Insert is used for the following products:

SDL Prod ID	Cat #	Description	
133	261	Shrink Seals (92X25mm)	400 seals/pkg
135	262	Shrink Seals (92X25mm)	6,000 seals/pkg
136	263	Shrink Seals - SPECIAL ORDER (92X25mm)	per gross (144)
138	265	Shrink Seals (special size)	
139	266	Shrink Seals-60mmX 25mm	1,000 seals/pkg
140	267	Shrink Seals 16 X 15	1,000 seals/pkg
141	268	Shrink Seals 18 X 15	1,000 seals/pkg
142	269	Shrink Seals 21 X 20	1,000 seals/pkg
143	270	Shrink Seals 27 X 25	1,000 seals/pkg
144	271	Shrink Seals 35 X 25	1,000 seals/pkg
145	272	Shrink Seals (BD Whiskey)	1,000 seals/pkg
146	273	Shrink Seals (BD 8A)	1,000 seals/pkg
147	274	Shrink Seals (ANAEROBIC)	1,000 seals/pkg
148	275	Shrink Seals (DIFCO EZDRAW)	1,000 seals/pkg
149	276	Shrink Seals (DIFCO ESP)	1,000 seals/pkg
150	277	Shrink Seals (ORGANON)	1,000 seals/pkg

Revision History

CR NUMBER	REVISION
0410-002	00
0910-002	01
1017-003	02
0420-001	03